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<http://www.servicebio.cn>

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Add: 5th Floor, C22 Building, Biopark, No. 388 Gaoxin 2nd Road, East Lake High-tech Developing Zone, Wuhan, Hubei, China, 430079 Website: <http://www.servicebio.com> Tel: 027-5111-3188

EMSA protocol

1. Experiment instruments

Name	Company	Model number
Microplate Reader	Rayto	RT-6100
Homogenizer	Servicebio	KZ-II
Centrifuge	DRAGONLAB	D3024R
Mini Centrifuge	Servicebio	D1008E
Vortex mixer	Servicebio	MX-F
Magnetic stirrers	Servicebio	MS-PB
Decolorizing shaker	Servicebio	TSY-B
Vertical Electrophoresis Apparatus	Servicebio	BV-2
Transfer Electrophoresis Apparatus	Servicebio	BT-2
Camera obscura	Servicebio	WGA0017
Bechtol	AIRTECH	SW-CJ-1FD
Darkroom lamp	Servicebio	WGA0016
Photographic film	Servicebio	WGJP0001
Image-analysis software	Adobe	Adobe PhotoShop
Gray analysis software	Alpha Innotech	alphaEaseFC
luminometer	CLINX	6300

2. Reagents

Reagents	Company	Catlog
Nuclear protein and Cytoplasmic Protein Extraction Kit	Beyotime	P0028
EMSA kits	Thermo Fisher	20148
10*TBE	Servicebio	G3002
40% Acrylamide	Servicebio	G2005
50% Glycerin	Servicebio	G0120
10% AP	Servicebio	G2003-5

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<http://www.servicebio.cn>

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Add: 5th Floor, C22 Building, Biopark, No. 388 Gaoxin 2nd Road, East Lake High-tech Developing Zone, Wuhan, Hubei, China, 430079 Website: <http://www.servicebio.com> Tel: 027-5111-3188

3. Protocol For EMSA

(1) **Histone extraction protocol for western blot** (Use Nuclear protein and cytoplasmic protein extraction kit)

- 1.1 **Prepare the reagents:** Dissolve the reagents at room temperature, then mix them up and put them on ice. The Nuclear protein and cytoplasmic protein extraction reagent A and B should be mixed up with PMSF until its final concentration comes to 1mM.
- 1.2 **Adherent cells:** Wash the cells three times with ice-cold PBS or EDTA. Centrifuge and collect the cells, discard the supernatant.
- 1.3 **Suspension cells:** Wash the cells three times with ice-cold PBS or EDTA. Centrifuge and collect the cells, discard the supernatant.
- 1.4 **From fresh tissues:** Cut the tissue into very small pieces as much as possible. Mix the appropriate amount of cytoplasmic protein extract A and B at a ratio of 20:1 (e.g., add 10 microliter of cytoplasmic protein extract B to 200 microliter of cytoplasmic protein extract A). Prepare a tissue homogenate with a final concentration of PMSF of 1 mM. Add 200μL of lysis buffer per 60mg tissue for machine homogenization. The homogenization should be performed in an ice bath or 4°C..Place in an ice bath for 15 minutes. After centrifugation at 2000g , 4°C for 5 minutes, the supernatant was transferred to a precooled EP tube, which was part of the cytoplasmic protein extracted.(Never touch the sediment when supernatant.)At this stage, some plasma proteins were still not extracted, and the remaining precipitation was again extracted by the method of cell precipitation.
- 1.5 Add 200μL of cytoplasmic protein extraction reagent A with PMSF into each 20μL of cell precipitate.
- 1.6 Vortex at the highest speed for 5s, let the cell precipitates completely suspended and dispersed;(Vortex time can be extended if cell precipitation were not completely suspended and dispersed).
- 1.7 Maintain constant agitation for 10-15 min in ice bath.
- 1.8 Add 10μL Cytoplasmic protein extraction reagent B(add PMSF before use).Vortex at the highest speed for 1s. Maintain constant agitation for 1 min in ice bath.
- 1.9 Vortex at the highest speed for 5s, Centrifuge at 12,000-16,000g, 4°C. for 5 min.
- 1.10 Aspirate the supernatant(plasma protein) and place in a fresh tube kept on ice.
- 1.11Aspirate the supernatant. Add 50μL nuclear protein extraction reagent B(add PMSF before use).
- 1.12 Vortex at the highest speed for 15-30s. Maintain constant agitation in ice bath. Vortex at

the highest speed for 15-30s in between 1-2 minutes.(Total:30min).

1.13 Centrifuge at 12,000-16,000g, 4°C for 10 min.

1.14 Transfer the supernatant to a fresh tube kept on ice. The supernatant is nuclear protein, you can immediately use or keep in -80°C fridge.

(2) Prepare 1.5mm native polyacrylamide gel and 0.5x TBE

Reagents	Concentration
	5.5%
H ₂ O ml	7.50
10*TBE μ l	500
40% Acrylamide ml	1.50
50% Glycerin μ l	500
10% AP μ l	50
TEMED μ l	10
Total ml	10

Pre-electrophoresis for 60 min at 120 V .Using precool 0.5×TBE as buffer solution.

(3) Prepare reaction mixtures as follows during pre-electrophoresis

Reagents	Negative Control μ l	Sample μ l	Competitor μ l	Mutated competitor μ l	Super Shift μ l
10×binding buffer	2	2	2	2	2
1 μ g/ μ l poly (dl.dc)	1	1	1	1	1
50% glycerol	1	1	1	1	1
1%NP-40	1	1	1	1	1
100mM MgCl ₂	1	1	1	1	1
200mM EDTA	1	1	1	1	1
Protein	0	2	2	2	2
Labeled Probe	2	2	2	2	2
Competitor	0	0	2	0	0
Mutated competitor	0	0	0	2	0
antibody	0	0	0	0	2
DEPC	Add to 20 μ l				

Place at room temperature for 25-30 min.

(4) Load the sample

After pre-electrophoresis ,change precool running buffer. Add 4 μ l of 6X DNA loading buffer to sample mixture. Load the samples on gel and run the gel at 150 V for 60min.

(5) Transfer

Put the positively charged nylon membrane into 0.5×TBE to equilibrate for 10 min. After the electrophoresis, remove the whole gel with the sample. The transferring condition is 300mA for 30min, using precool 0.5×TBE as buffer solution

(6) Crosslinking

Place the membrane containing the reaction mixture under the UV lamp (20 cm below the light) for 20min.

(7) Blocking

Block the membrane on a table concentrator for 20min at room temperature using blocking buffer.

(8) Antigen-Antibody Reaction

Aspirate the blocking buffer. Incubate the membrane with appropriate dilutions of primary antibody (1:300). for 30min at room temperature.

(9) Washing

Wash the membrane for 4 times(10min each) with 1X washing buffer.

(10) Chemiluminescence:

Take the membrane out of the eluent, absorb the liquid, put it in a chemiluminescence instrument, add ECL, wait for 1-2 minutes, absorb the excess liquid, and start chemiluminescence.

4. Result and analysis (shown in Excel)