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2. Reagents

Reagents	Company	Catlog
IP Lysis Buffer	Servicebio	G2038-100ML
Protein A/G-beads	Millipore	IP05
50*cocktail	Servicebio	G2006
PMSF (100mM)	Servicebio	G2008
Phosphatase Inhibitor	Servicebio	G2007
BCA protein quantitative detection kit	Servicebio	G2026
5*/2*loading buffer	Servicebio	G2013/G2031
SDS-PAGE kits	Servicebio	G2003
Protein Marker	Servicebio	26617

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Add: 5th Floor, C22 Building, Biopark, No. 388 Gaoxin 2nd Road, East Lake High-tech Developing Zone, Wuhan, Hubei, China, 430079 Website: <http://www.servicebio.com> Tel: 027-5111-3188

PVDF membrane 0.45μm	Servicebio	G6015-0.45
Milk	Servicebio	G5002-100G
TWEEN 20	Servicebio	WGT8220
ECL	Servicebio	G2014
Developer and fixer	Servicebio	G2019
PeroXidase-conjugated Goat Anti-Rabbit IgG(H+L)	Servicebio	GB23303
PeroXidase-conjugated Donkey Anti-Goat IgG(H+L)	Servicebio	GB23404
PeroXidase-conjugated Goat Anti-Mouse IgG(H+L)	Servicebio	GB23301
PeroXidase-conjugated Goat Anti-Rat IgG(H+L)	Servicebio	GB23302
Transfer buffer	Servicebio	G2017
Electrophoretic Buffer	Servicebio	G2018
TBS buffer solution	Servicebio	G0001

3. Protocol For CoIP

(1) Cell protein extraction

1.1 Cell culture (Adherent cells: The cells should be about 80% - 90% confluent in a 100 mm cell culture dish, or the number of cells in the culture bottle was $2-5 \times 10^7$. Suspension cells: Centrifuge to collect cells ($2-5 \times 10^7$))

1.2 Add enough cold IP buffer (Protease and phosphatase inhibitors should be added before use). Maintain constant agitation for 10 min in ice bath, use pipette to mix it up. Then gently transfer the cell suspension into a 1.5ml microcentrifuge tube. Continue splitting for 20 min in ice bath.

1.3 Centrifuge at 12,000 rpm, 4°C for 10 min. Gently remove the tubes from the centrifuge and place on ice, aspirate the supernatant and place in a fresh tube kept on ice; then determination with the BCA kit

1.4 Denature some lysis solution for Input

(2) Tissue protein extraction

2.1 Wash the sample with cold PBS for twice or three times, then homogenize with an electric homogenizer. Volumes of lysis buffer must be determined in relation to the amount of tissue present.

2.2 Gently transfer the cell suspension into a 1.5ml microcentrifuge tube. Maintain constant agitation for 30 min in ice bath, use pipette to mix it up.

2.3 Centrifuge at 12,000 rpm, 4°C for 10 min. Gently remove the tubes from the centrifuge and place on ice, aspirate the supernatant and place in a fresh tube kept on ice; Use BCA method to determine protein concentration.

2.4 Denature some lysis solution for Input.

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Wuhan Servicebio Technology Co., Ltd.

3.1 Add 1.0 μ g IgG and 20 μ L Protein A/G PLUS-agarose to the IgG group. Add 20 μ L Protein A/G PLUS-agarose to the experimental group. Shaken at 4 °C for 1 h.

3.2 Centrifuge at 2000g, 4°C for 5 min. Gently remove the tubes from the centrifuge and place on ice, aspirate the supernatant.

3.3 Incubate 1-10 μ L (0.2-2 μ g) antibody coupled beads in 250 μ L cell lysis supernatants, rotate overnight at 4 °C

3.4 Mix the lysates with 80 μ L Protein A/G PLUS-agarose rotating at 4 °C for 2 h.

3.5 Centrifuge at 2000g, 4°C for 5 min. gently remove the supernatant.

3.6 Wash the immunoprecipitation complex with 1 ml of precooled IP lysate (without any inhibitors) for 4 times. Centrifuge at 2000g, 4°C for 10 min, and the supernatant should be carefully discarded each time.

3.7 After the last washing, discard the supernatant as much as possible, then add 80 μ L 1 $\times$ reduced loading buffer, boil for 10 min in boiling water. Centrifuge at 1000 g, 4 °C, for 5 min. then aspirate the supernatant.

4. Detection

4.1 Separate the protein sample by SDS-PAGE (5-15%) and analyze by western blot.

4.2 The results were analyzed in the form of Excel.