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Paraffin-fluorescence Probe-FISH and Immunofluorescence Protocol

1. Apparatus and Reagents

1.1 Apparatus

Name	Producer	Model
Dehydrator	WHJJ	JJ-12J
Paraffin embedding machine	WHJJ	JB-P5
Pathologic microtome	Leica	RM2016
Frozen flat	WHJJ	JB-L5
Water Bath-Slide	Kedee	KD-P
RNase-free glass microscope slides	Servicebio	
Micro-centrifuge	Servicebio	D1008E
Rocker	Servicebio	TSY-B
Vortex	Servicebio	MV-100
Pipettor	Dragon	KE0003087/KA0056573
Microscopy	Nikon	NIKON ECLIPSE CI
Imaging system	Nikon	NIKON DS-U3
Liquid Blocker PAP Pen	Servicebio	G6100
Refrigerator	HAIER	BCD-192TGN
Incubator	LABOTERY	GSP-70
Autoclav	PANASONIC	MLS-3751L-PC

1.2 Major reagents

Reagent	Manufacturer	Article Number
Ethanol	SCRC	100092683
BioDewax and clear solution	Servicebio	G1128
PBS	Servicebio	G0002
Proteinase K	Servicebio	G1234
4% of paraformaldehyde (DEPC water)	Servicebio	G1113
DEPC	Amresco	E174
20×SSC solution	Servicebio	G3016-4

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Add: 5th Floor, C22 Building, Biopark, No. 388 Gaoxin 2nd Road, East Lake High-tech Developing Zone, Wuhan, Hubei, China, 430079 Website: <http://www.servicebio.com> Tel: 027-5111-3188

BSA	Servicebio	G5001
DAPI	Servicebio	G1012
Anti-fluorescence quenching	Servicebio	G1401
sealing tablets		
Hybridization buffer	Servicebio	G3016-3
First antibody		
Second antibody		

2. The steps of the experiment

2.1 **Organization fixation:** take out the organization, wash clean, then Immediately put in the fixed fluid (DEPC) above 12h.

2.2 **Dehydration:** the tissue is dehydrated by gradient alcohol, paraffin, embedding.

2.3 **Section:** the paraffin is sliced through the slicer, the piece of the slice machine and the 62 - degree oven roast for 2 hours.

2.4 **Dewaxing and dehydration:** soak sections in 2 changes of Dewaxing Transparent Liquid, 15 minutes each. Dehydrate in 2 changes of pure ethanol for 5 minutes each. Then, followed respectively by dehydrating in gradient ethanol of 85%and 75% ethanol 5 minutes each. Wash in DEPC dilution.

2.5 **Digestion:** according to the tissue fixation time, the slices are boiled in the retrieval solution for 10-15 minutes and naturally cooled. Mark the objective tissue with liquid blocker pen, according to the characteristics of tissues, Add proteinase K(20ug/ml) working solution to cover objectives and incubate at 37°C for___min. Wash in pure water, then wash three times in PBS (pH 7.4) on a Rocker device, 5 min each.

2.6 **Pre-hybridization:** add Pre-hybridization solution to each section and incubate for 1 h at 37°C.

2.7 **Hybridization:** remove the pre-hybridization solution, add the probe hybridization solution with concentration of ____, and incubate the section in a humidity chamber and hybridize overnight at___°C.

2.8 **Washing:** remove the hybridization solution. Wash sections in 2×SSC for 10 min at 37°C , wash sections in 1×SSC two times for 5 min each at 37°C, and wash in 0.5×SSC for 10 min at room temperature. Formamide washing can be added if there are more non-specific hybrids.

2.9 **Blocking:** add blocking serum to the section and incubate at room temperature for 30 min.

2.10 **Incubate first antibody:** PBS solution containing a____dilution of primary antibody were added and incubated at 4°C overnight. Samples were then washed with PBS three times for 5 min each at RT.

2.11 **Incubate second antibody:** after washing, the section was incubated for 50 min with second antibody at RT. Samples were then washed with PBS three times for 5 min each at RT.

2.12 **Stain cell nuclei (counter stain):** incubate with DAPI for 8min in the dark, and then mounting.

2.13 **Microscopic examination and photography:** to take photos with positive fluorescence microscope. DAPI glows blue by UV excitation wavelength 330-380 nm and emission wavelength 420 nm; FAM glows green by excitation wavelength 465-495 nm and emission wavelength 515-555 nm; CY3 glows red by excitation wavelength 510-560 nm and emission wavelength 590 nm.

3. Interpretation of the results

The nuclear stained by DAPI were blue under ultraviolet excitation, and the positive expression was a kind of fluorescence labeled by corresponding luciferin. FAM (488) is green light, cy3 is red light. The results of mRNA in situ hybridization were cytoplasmic positive and a few nuclear positive were normal. MicRNA and lncRNA were expressed differently. According to the expression, different fluorescence brightness is strong or weak.

Attached table 1 probe information.

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Wuhan Servicebio Technology Co., Ltd.