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Add: 5th Floor, C22 Building, Biopark, No. 388 Gaoxin 2nd Road, East Lake High-tech Developing Zone, Wuhan, Hubei, China, 430079 Website: <http://www.servicebio.com> Tel: 027-5111-3188

BSA	Servicebio	G5001
蛋白酶 K	Servicebio	G1234
核固红染液	Servicebio	G1035
甘油明胶	Servicebio	G1402
杂交缓冲液	Servicebio	G3016-3
anti-DIG-AP	jackson	200-052-156
BCIP/NBT 显色剂	博士德	AR1023

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1.2 Major reagents

Reagent	Manufacturer	Article Number
Ethanol	SCRC	100092683
BioDewax and Clear solution	Servicebio	G1128
PBS	Servicebio	G0002
Proteinase K	Servicebio	G1234
4% of paraformaldehyde (DEPC water)	Servicebio	G1113
DEPC	Amresco	E174
20×SSC solution	Servicebio	G3016-4
BSA	Servicebio	G5001
Nuclear fast red dye	Servicebio	G1035
Glycerin-gelatin	Servicebio	G1402
hybridization buffer	Servicebio	G3016-3

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Add: 5th Floor, C22 Building, Biopark, No. 388 Gaoxin 2nd Road, East Lake High-tech Developing Zone,
Wuhan, Hubei, China, 430079 Website: <http://www.servicebio.com> Tel: 027-5111-3188

anti-DIG-AP	Jackson	200-052-156
NBT/BCIP reagent	Boster	AR1023

2. The Steps of the Experiment

2.1. **Organization fixed:** Take out the organization, wash clean, then Immediately put in the fixed fluid (DEPC) to fix 2-12h.

2.2. **Dehydration:** The tissue is dehydrated by gradient alcohol and then paraffin embedded.

2.3. **Section:** The paraffin is sliced through the slicer, the piece of the slice machine and the 62 - degree oven roast for 2 hours.

2.4. **Dewaxing and dehydration:** Soak sections in 2 changes of Dewaxing Transparent Liquid, 15 minutes each. Dehydrate in 2 changes of pure ethanol for 5 minutes each. Then, followed respectively by dehydrating in gradient ethanol of 85%and 75% ethanol 5 minutes each. Wash in DEPC dilution.

2.5. **Digestion:** According to the tissue fixation time, the slices are boiled in the retrieval solution for 10-15 minutes and naturally cooled. Mark the objective tissue with liquid blocker pen, according to the characteristics of tissues, Add proteinase K(20ug/ml) working solution to cover objectives and incubate at 37°C for___min. Washing with pure water, then wash three times with PBS (pH 7.4) in a Rocker device, 5 min each.

2.6. **Pre-hybridization:** Add Pre-hybridization solution to each section and incubate for 1 h at 37°C.

2.7. **Hybridization:** Discard the pre-hybridization solution, add the probe hybridization solution ,concentration____, and incubate the section in a humidity chamber and hybridize overnight at___°C.

2.8. **Washing:** Remove the hybridization solution. Wash sections with 2×SSC for 5 min at 37°C, Wash sections in 1×SSC two times for 5 min each at 37°C, and wash in 0.5×SSC for 10 min at room temperature. Formamide washing can be added if there are more non-specific hybrids.

2.9. **Blocking:** add blocking solution (Rabbit serum) to the section and incubate at room temperature for 30 min.

2.10. Add the Alkaline Phosphatase-conjugated IgG Fraction Monoclonal Mouse Anti-Digoxin Antibody: (anti-DIG-AP): Remove the blocking solution and add anti-DIG-HRP. Incubate at 37 °C for 40 min, then Wash sections in TBS four times for 5 min each.

2.11. **NBT/BCIP developing:** dry sections slightly, add fresh prepared NBT/BCIP chromogenic reagent to marked tissue. Manage reaction time by observing in microscopy until positive expression appears shows blue-purple. Then stop developing reaction by wash in running tap water.

- 2.12. **Multicolor nucleus:** dripping with red dye, and then rinse with pure water.
- 2.13. **Mounting:** Natural air dry and mount with Glycerol Jelly mounting medium.
- 2.14. **Image the results using a brightfield fluorescence microscope.** Store slides at RT.

3. Interpretation of the Results

The nuclei were red, and the positive expression of BCIP/NBT was blue-purple. The results of mRNA in situ hybridization showed that cytoplasm was positive and a few nuclear positive were normal. The expression of micRNA and lncRNA was different. According to the expression, the color is very light.

Attached table 1 probe information.

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Wuhan Servicebio Technology Co., Ltd.