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Frozen Tissue-fluorescence Probe-FISH Protocol

1. Apparatus and Reagents

1.1 Apparatus

Name	Producer	Model
Frozen slicer	Thermo	Cryotome E
Glass slide	Servicebio	
Coverslip	Citotest	10212432C
Enzyme-free centrifuge tubes	Servicebio	EP-150-M
Shaker	Servicebio	TSY-B
Vortex	Servicebio	MV-100
Pipettor	Dragon	KE0003087/KA0056573
Liquid Blocker PAP Pen	Servicebio	G6100
Refrigerator	HAIER	BCD-192TGN
Microscopy	Nikon	NIKON ECLIPSE CI
Imaging system	Nikon	NIKON DS-U3
Incubator	LABOTERY	GSP-70
autoclave	PANASONIC	MLS-3751L-PC

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Add: 5th Floor, C22 Building, Biopark, No. 388 Gaoxin 2nd Road, East Lake High-tech Developing Zone, Wuhan, Hubei, China, 430079 Website: <http://www.servicebio.com> Tel: 027-5111-3188

2. The Steps of the Experiment

- 2.1. **Tissue fixation:** tissue was washed and fixed with 4% paraformaldehyde for 12 h or more.
- 2.2. **Tissue dehydration:** Transfer tissue to 15% sucrose solution for 8 h and were subsequently to a 30% sucrose solution overnight.
- 2.3. **Frozen section fixation:** The frozen sections were dried at room temperature, fixed in 4% paraformaldehyde (DEPC) for 10 min, and shaken on the shaker in PBS (pH 7.4) three times for 5 min each.
- 2.4. **Digestion:** Mark the objective tissue with liquid blocker pen, according to the characteristics of different tissues and different indicators, add proteinase K (20 ug/ml) to cover tissues and incubate at 37°C for ____ minutes. Washed with sterilized water and then washed three times with PBS for 5 min each time.
- 2.5. **Prehybridization:** add hybridization buffer onto specimen and incubate at 37°C for 1h.
- 2.6. **Hybridization:** remove the pre-hybridization solution, add the probe hybridization solution with concentration of ____, and incubate the section in a humidity chamber and hybridize overnight at ____°C.
- 2.7. **Washing:** Remove the hybridization solution. Wash sections in 2×SSC for 10 min at 37°C , 1×SSC two times for 5 min each at 37°C, and 0.5×SSC for 10 min at room temperature. Formamide washing can be added if there are more non-specific hybrids.
- 2.8. **Stain cell nuclei (counter stain):** Incubate with DAPI for 8min in the dark, and then mounting with anti-fluorescence quenching sealing tablets.
- 2.9. **Microscopic examination and photography :** To take photos with positive fluorescence microscope. DAPI glows blue by UV excitation wavelength 330-380 nm and emission wavelength 420 nm; FAM glows green by excitation wavelength 465-495 nm and emission wavelength 515-555 nm; CY3 glows red by excitation wavelength 510-560 nm and emission wavelength 590 nm.

3. Interpretation of the Results

The nuclear stained by DAPI were blue under ultraviolet excitation, and the positive expression was a kind of fluorescence labeled by corresponding luciferin. FAM (488) appears green, cy3 appears red. The results of mRNA in situ hybridization were cytoplasmic positive and a few nuclear positive were normal. MicRNA and lncRNA were expressed differently. According to the expression, Different fluorescence brightness is strong or weak.

Note: All reagents, instrument need RNase free.

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Wuhan Servicebio Technology Co., Ltd.