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Add: 5th Floor, C22 Building, Biopark, No. 388 Gaoxin 2nd Road, East Lake High-tech Developing Zone,
Wuhan, Hubei, China, 430079 Website: <http://www.servicebio.com> Tel: 027-5111-3188

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NBT/BCIP-DIG-ISH Protocol (Frozen Tissue)

1. Apparatus and Reagents

1.1 Apparatus

Name	Producer	Model
Frozen slicer	Thermo	Cryotome E
Glass slide	Servicebio	
Coverslip	Citotest	10212432C
Enzyme-free centrifuge tubes	Servicebio	EP-150-M
Shaker	Servicebio	TSY-B
Vortex	Servicebio	MV-100
Pipettor	Dragon	KE0003087/KA0056573
Liquid Blocker PAP Pen	Servicebio	G6100
Refrigerator	HAIER	BCD-192TGN
Microscopy	Nikon	NIKON ECLIPSE CI
Imaging system	Nikon	NIKON DS-U3
Incubator	LABOTERY	GSP-70
autoclave	PANASONIC	MLS-3751L-PC

1.2 Major reagents

Reagent	Manufacturer	Article Number
DEPC	Amresco	E174
4% of paraformaldehyde (DEPC water)	Servicebio	G1113
OCT embedding agent	Servicebio	G6059-110ML
Ethanol	SCRC	100092683
Xylene	SCRC	10023418
Sucrose	Servicebio	G5031
PBS	Servicebio	G0002
20×SSC solution	Servicebio	G3016-4
BSA	Servicebio	G5001
Proteinase K	Servicebio	G1234
Nuclear fast red dye	Servicebio	G1035

Glycerin-gelatin	Servicebio	G1402
Hybridization buffer	Servicebio	G3016-3
anti-DIG-AP	Jackson	200-052-156
NBT/BCIP reagent	Boster	AR1023

2. The Steps of the Experiment

2.1. **Tissue fixation:** take out the organization, wash clean, then Immediately put in the fixed fluid (DEPC) above 12h.

2.2. **Tissue dehydration:** Transfer tissue to 15% sucrose solution for 8 h and were subsequently to a 30% sucrose solution overnight.

2.3. **Frozen section fixation:** The frozen sections were dried at room temperature, fixed in 4% paraformaldehyde (DEPC) for 10 min, and shaken on the shaker in PBS (pH 7.4) for three times, 5 min each.

2.4. **Digestion:** Mark the objective tissue with liquid blocker pen, according to the characteristics of different tissues and different indicators, add proteinase K (20 ug/ml) to cover tissues and incubate at 37°C for _____minutes. Washed with sterilized water and then washed three times in PBS for 5 min each time.

2.5. **Prehybridization:** add hybridization buffer onto specimen and incubate at 37°C for 1h.

2.6. **Hybridization:** remove the pre-hybridization solution, add the _____ probe hybridization solution with concentration of _____, and incubate the section in a humidity chamber and hybridize overnight at _____°C.

2.7. **Washing:** Remove the hybridization solution. Wash sections in 2×SSC for 10 min at 37°C, 1×SSC two times for 5 min each at 37°C, and 0.5×SSC for 10 min at room temperature. Formamide washing can be added if there are more non-specific hybrids.

2.8. **Blocking:** add blocking solution (Rabbit serum) to the section and incubate at room temperature for 30 min.

2.9. **Add the Alkaline Phosphatase-conjugated IgG Fraction Monoclonal Mouse Anti-Digoxin Antibody:** (anti-DIG-AP): Remove the blocking solution and add anti-DIG-HRP. Incubate at 37 °C for 40 min, then Wash sections in TBS four times for 5 min each.

2.10. **NBT/BCIP developing:** dry sections slightly, add fresh prepared NBT/BCIP chromogenic reagent to marked tissue. Manage reaction time by observing in microscopy until positive expression appears shows blue-purple. Then stop developing reaction by wash in running tap water.

2.11. **Multicolor nucleus:** dripping with red dye, and then rinse with pure water.

2.12. **Mounting:** Dehydrate successively in gradient ethanol of 75%, 85%, and 2 changes of pure

ethanol, respectively, 6 min each. Clear in xylene for 6 min and mount with resin mounting medium.

2.13. Image the results using a brightfield fluorescence microscope. Store slides at RT.

3. Interpretation of the Results

The nucleus stained with hematoxylin are blue, and the positive expression of DAB was brownish yellow. The results of in situ hybridization of mRNA were cytoplasmic positive and a few nuclear positive were normal. MicRNA and lncRNA were expressed differently. According to the expression of color depth.

Note: All reagents, instrument need RNase free.

Attached table 1 probe information.