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Cell Climbing Slides EDU Staining Protocol

1. Apparatus and Reagents

1.1 Apparatus

Name	Producer	Model
Glass microscope slides	Servicebio	
Vortex	Servicebio	MX-F
Micro-centrifuge	Servicebio	D1008E
Rocker	Servicebio	TSY-B
Pipettor	Dragon	KE0003087/KA0056573
Liquid blocker pen	Servicebio	G6100
Ortho-Fluorescent Microscopy	Nikon	NIKON ECLIPSE C1
Imaging system	Nikon	NIKON DS-U3

1.2 Major reagents

Name	Producer	Code	Dilution
PBS solution	Servicebio	G0002	
Permeabilize solution	Servicebio	G1204	
Glycine	Servicebio	G5010-500G	
EDU staining kit	Servicebio	G1601	
		G1602	
		G1603	
		G1604	
DAPI	Servicebio		
Anti-fade mounting medium	Servicebio	G1401	

2. Procedure

2.1 Add 100 μ L of penetrant to each well (0.5% TritonX-100 of PBS), incubate in rocker device for 10 minutes; wash with PBS once for 5 minutes.

2.2 Add 100 μ L EDU staining reaction solution to each well, and discard the staining reaction solution after incubating for 30 minutes at room temperature and protected from light;

Component	Volume
Reaction buffer	935ul
Catalytic reagent (A)	10ul
Fluorescent dye(B)	5ul
Catalytic additive(C)	50ul
Total volume	1000ul

2.3 DAPI counterstain in nucleus: Wash three times with PBS (pH 7.4) in a Rocker device, 5 min each. Then incubate with DAPI solution at room temperature for 10 min, kept in dark place.

2.4 Mount: Throw away liquid slightly, then coverslip with anti-fade mounting medium

2.5 Microscopy detection and collect images by Fluorescent Microscopy. DAPI glows blue by UV excitation wavelength 330-380 nm and emission wavelength 420 nm; FITC glows green by excitation wavelength 465-495 nm and emission wavelength 515-555 nm; 555 glows red by excitation wavelength 550 nm and emission wavelength 590 nm.

3. Results

Nucleus is blue by labeling with DAPI. Positive cells are green or red according to the fluorescent labels used.