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Add: 5th Floor, C22 Building, Biopark, No. 388 Gaoxin 2nd Road, East Lake High-tech Developing Zone, Wuhan, Hubei, China, 430079 Website: <http://www.servicebio.com> Tel: 027-5111-3188

Phalloidine Staining Protocol for Cell Sample (Red)

1. Apparatus and Reagents

1.1 Apparatus

Name	Producer	Model
Glass microscope slides	Servicebio	
Vortex	Servicebio	MX-F
Micro-centrifuge	Servicebio	D1008E
Rocker	Servicebio	TSY-B
Pipettor	Dragon	KE0003087/KA0056573
Liquid blocker pen	Servicebio	G6100
Ortho-Fluorescent Microscopy	Nikon	NIKON ECLIPSE C1
Imaging system	Nikon	NIKON DS-U3

1.2 Major reagents

Name	Producer	Code	Dilution
PBS solution	Servicebio	G0002	
Permeabilize solution	Servicebio	G1204	
phalloidin (Cy3)	Servicebio	G1041	1:200
DAPI	Servicebio	G1012	
anti-fade mounting medium	Servicebio	G1401	

2. Procedure

2.1 Permeabilization: Dry the cell climbing slides slightly. Mark the objective area with liquid blocker pen, where add 50-100 μ l of permeabilize working solution. Incubate for 20 min at room temperature. Wash three times with PBS solution, 5 min each.

2.2 Phalloidine (Cy3) stain: Add phalloidine diluted solution to cover with the objective area, incubate at room temperature for 2 h in the dark.

2.3 DAPI counterstain in nucleus: Wash three times with PBS (pH 7.4) in a Rocker device, 5 min each. Then incubate with DAPI solution at room temperature for 10 min in the dark.

2.3 Mount: Wash three times with PBS (pH 7.4) in a Rocker device, 5 min each. Throw away liquid slightly. Put the slides on a glass microscope slide and then mount with anti-fade mounting medium.

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Wuhan Servicebio Technology Co., Ltd.

2.4 Microscopy detection and collect images by Fluorescent Microscopy. DAPI glows blue by UV excitation wavelength 330-380 nm and emission wavelength 420 nm; FITC glows green by excitation wavelength 465-495 nm and emission wavelength 515-555 nm; CY3 glows red by excitation wavelength 510-560 nm and emission wavelength 590 nm.

3. Results

Nucleus is blue by labeling with DAPI. Positive cells labelled by Phalloidine (Cy3) are red.