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Add: 5th Floor, C22 Building, Biopark, No. 388 Gaoxin 2nd Road, East Lake High-tech Developing Zone, Wuhan, Hubei, China, 430079 Website: <http://www.servicebio.com> Tel: 027-5111-3188

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PASM Staining Experimental Report

1. Lab Equipment and Reagents

1.1 Lab equipment

Items	Manufacturer	Model
Dehydrator	DIAPATH	Donatello
embedding machine	Wuhan Junjie Electronics Co., Ltd.	JB-P5
Pathology microtome	Shanghai Leica Instruments Co., Ltd.	RM2016
Frozen platform	Wuhan Junjie Electronics Co., Ltd.	JB-L5
Water Bath-Slide Drier	Zhejiang Jinhua Kedi Instrumental Equipment CO.,LTD	KD-P
Laboratory oven	Tianjin Labotery Instrument Equipment Co., Ltd.	GFL-230
Microscope slide	Servicebio	
Upright optical microscope	Nikon Japan	Nikon Eclipse E100
Imaging system	Nikon Japan	NIKON DS-U3

1.2 Chemical reagents

Items	Manufacturer	Model
Absolute alcohol	Sinopharm Chemical Reagent Co., Ltd.	100092683
Xylene	Sinopharm Chemical Reagent Co., Ltd.	10023418
PASM staining kit	Servicebio	G1059
Neutral balsam	Sinopharm Chemical Reagent Co., Ltd.	10004160

2. Experimental Steps

2.1 Paraffin section deparaffinization and rehydration: put the slides into xylene I 20minutes-xylene II 20 minutes-absolute ethanol I 5 min-absolute ethanol II 5 min-75% alcohol for 5 min, tap water washing and then distilled water washing.

2.2 Periodic acid-silver metharamine staining: slides into periodic acid-silver metharamine (PASM) staining solution A, soak overnight, tap water washing.

2.3 Working solution preparation: PASM staining solution D: PASM staining solution C 20:1 mixed into PASM stock solution, take PASM stock solution 20ml, adding 15ml of distilled water,

And then add PASM staining solution E 2ml and well mixed.

2.4 Tissue acidification: put the slides into PAS staining solution B for 15 to 20 minutes, and distilled water washing for 3-5 times;

2.5 Incubation and staining: the working solution placed in a 60°C oven for preheating in advance, and add the preheated working solution in the slides drop by drop, incubation with cover (56°C-59°C);

2.6 Color rendering: after 40 minutes of incubation, take out the slides and observe under a microscope (the working solution doesn't need to be washed away during observation), and the kidney is subject to the glomerular capillary basement membrane being black. If the coloring is not deep enough, continue to incubate in the oven; Lungs are subject to fungus turning black (if there is no fungus, the bronchial darkening is the standard). After that, take out the slides and observe under a microscope every 10 minutes until staining result is okay, distilled water washing for three times.

2.7 PASM staining solution F treatment: after treatment of PASM staining solution F, wash once with distilled water.

2.8 Counterstaining PASM staining solution F (optional): 85%, 95% alcohol for gradient dehydration, put the slides into PASM staining solution F for 20 seconds.

2.9 Dehydration and sealing: put the slides into absolute ethanol I 5min- absolute ethanol II 5min- absolute ethanol III 5min- xylene I 5min-xylene II 5min transparent, neutral balsam sealing.

2.10 Microscope examination, image collection and analysis.

3. Results

The glomerular basement membrane, the renal capillary basement membrane and the renal tubular epithelial basement membrane are black, and the fungi and collagen fibers are also grayish black or brownish black with red background.

4. Note

4.1 The slides should be cleaned with distilled water before incubation, and the container for preparing the PASM stock solution and working solution must be cleaned, otherwise, impurities will easily appear.

4.2 PASM stock solution and working solution should be ready for use.