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Add: 5th Floor, C22 Building, Biopark, No. 388 Gaoxin 2nd Road, East Lake High-tech Developing Zone, Wuhan, Hubei, China, 430079 Website: <http://www.servicebio.com> Tel: 027-5111-3188

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## Safranin Fast Green Staining for Bone Experimental Report

### 1. Experimental Equipment and Reagents

#### 1.1 Experimental equipment

| Name                         | Manufactor                                       | Model              |
|------------------------------|--------------------------------------------------|--------------------|
| Dehydrator                   | DIAPATH                                          | Donatello          |
| Embedding machine            | Wuhan JunjieElectronics Co., Ltd                 | JB-P5              |
| Pathological section machine | Shanghai Leica Instrument Co., Ltd               | RM2016             |
| Frozen platform              | Wuhan Junjie Electronics Co., Ltd                | JB-L5              |
| KD-P Water Bath              | Kedee                                            | KD-P               |
| Oven                         | Tianjin Lai Bo Rui Instrument Equipment Co., Ltd | GFL-230            |
| Slides                       | Servicebio                                       |                    |
| Orthostatic microscope       | Nikon, JAPAN                                     | Nikon Eclipse E100 |
| Image system                 | NIKON, JAPAN                                     | Nikon DS-U3        |

#### 1.2 Main experimental reagents

| Reagent name                | Manufactor                               | Article number |
|-----------------------------|------------------------------------------|----------------|
| Absolute ethanol            | Sinopharm Group Chemical Reagent Co. LTD | 100092683      |
| Xylene                      | Sinopharm Group Chemical Reagent Co. LTD | 10023418       |
|                             | Servicebio                               |                |
| BioDewax and Clear Solution |                                          | G1128          |
| SaffronFastGreen dye        | Servicebio                               | G1053          |
| solution kit(Bone Tissue)   |                                          |                |
| Differentiation solution    | Servicebio                               | G1040          |
| Neutral resin               | Sinopharm Group Chemical Reagent Co. LTD | 10004160       |

### 2. Experimental Steps

**2.1** Paraffin slides dewaxed as follow: Two changes of pure BioDewax and Clear Solution for 20 min. Two changes of pure ethanol for 5min. 75% ethanol for 5min. Keep slides in tap water.

**2.2** Fast green staining: The slides stained in fast green dye solution for 1-5 minutes, washed away the excess dye solution until the cartilage was colorless, and soaked in 1% hydrochloric acid and alcohol for 10s. And then washing them with tap water.

**2.3 Saffron staining:** The slides were stained in saffron dye solution for 1-5s, and then put into four cylinders of absolute ethanol, where three for rapid dehydration for 5s, 2s, 10s respectively, and kept in the fourth cylinder.

**2.4** The slides were immersed in xylene to transparent for 5min, sealing with neutral resin.

**2.5** Observed under microscope, and took images for analysis.

### **3. The Results Were As Follows**

Cartilage was red or orange-red, and bone formation was green. And some connective tissues were red.

### **4. Precautions**

4.1 During the dyeing process, the saffron staining cannot be overtime, because it is easy to hybridize with fast green to become purple-blue.

4.2 If the result of the staining is not well, it can be washed with hydrochloric acid solution and alcohol until the color fades, and then re-stained.

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Wuhan Servicebio Technology Co., Ltd.