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Bone Tissue Goldner Trichrome Staining

Experiment Report

1. Experimental Equipment and Reagents

1.1 Experimental equipment

Equipment name	Manufacturer	Model No.
Dehydrator	DIAPATH	Donatello
Embedding center	Wuhan Junjie Electronics Co., Ltd.	JB-P5
Pathological microtome	Shanghai Leica Instrument Co., Ltd.	RM2016
Cooling plate	Wuhan Junjie Electronics Co., Ltd.	JB-L5
Tissue spreading water bath	Zhejiang Jinhua Kedi Instrumental Equipment Co., Ltd.	KD-P
Oven	Tianjin Leibo Terry Equipment Co., Ltd.	GFL-230
Microscope slide	Servicebio	
Upright electron microscope	JAPAN NIKON	NIKON ECLIPSE E100
Imaging system	JAPAN NIKON	NIKON DS-U3

1.2 Main experiment reagents

Reagent name	Manufacturer	Item No.
Anhydrous ethanol	Sinopharm Chemical Reagent Co.,Ltd.	100092683
Xylene	Sinopharm Chemical Reagent Co.,Ltd.	10023418
Goldner staining solution suit	Servicebio	G1064
Glacial acetic acid	Sinopharm Chemical Reagent Co.,Ltd.	10000218
Neutral balsam	Sinopharm Chemical Reagent Co.,Ltd.	10004160

2. Experimental Procedure

2.1 Deparaffinize the paraffin sections to water: Place the sections in sequence in Xylene I for 20min - Xylene II for 20min - Anhydrous ethanol I for 5min - Anhydrous ethanol II for 5min - 75% Ethyl alcohol for 5min , wash with tap water.

2.2 Nuclear staining: Goldner staining solution A and Goldner staining solution B are mixed in equal proportions, put sections in the mixed solution and stain for 20 min, wash with tap water, then differentiate with 1% hydrochloric acid alcohol solution for 2s, resin with tap water, and wash with distilled water.

2.3 Staining with Goldner staining solution C: Immerse sections into first jar of Goldner staining solution C to stain for 5-10 minutes, then rinse quickly with 0.2% glacial acetic acid, for 3 seconds per tank.

2.4 Staining with Goldner staining solution D: Dip the sections into Goldner staining solution D to stain for 3 min, control under microscope, subject to the fading of Goldner staining solution C at the collagen area.

2.5 Staining with Goldner staining solution C: Dip the sections in second jar of Goldner staining solution C and stain for 3-5min, rinse quickly with 0.2% glacial acetic acid for 3s per jar.

2.6 Staining with Goldner staining solution E: Put the sections into Goldner staining solution E and stain for 3-5min, the differentiate with three jars of 0.2% glacial acetic acid, for 2s per jar, and dehydrate with three jars of anhydrous ethanol, for 2s, 3s and 5s respectively.

2.7 Transparency and sealing: Place the sections into the third jar of anhydrous ethanol for 5min, Xylene transparency for 5min, lastly, seal with neutral balsam.

2.8 Microscope inspection, image acquisition and analysis.

3. Interpretation of Results:

The cell nucleus are blue-violet, the mineralized bone is green, and the non-mineralized bone is orange-red.

4. Precautions:

4.1 According to the number of sections, replace the iron hematoxylin staining solution regularly and it can also be prepared when using.

4.2 The staining degree of Goldner staining solution D and Goldner staining solution C must be well controlled. The collagen part cannot be too red, it will affect the coloring of Goldner staining solution D.

4.3 During the differentiating of Goldner staining solution E by glacial acetic acid, the differentiation should be proper, and not be too deep or too shallow.

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Wuhan Servicebio Technology Co., Ltd.