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Gram Staining Experiment Result

I. Experimental equipment and reagents

1. Experimental equipment

Equipment name	Manufacturer	Model No.
Dehydrator	DIAPATH	Donatello
Embedding center	Wuhan Junjie Electronics Co., Ltd.	JB-P5
Pathological microtome	Shanghai Leica Instrument Co., Ltd.	RM2016
Cooling plate	Wuhan Junjie Electronics Co., Ltd.	JB-L5
Tissue spreading water bath	Zhejiang Jinhua Kedi Instrumental Equipment Co., Ltd.	KD-P
Oven	Tianjin Leibo Terry Equipment Co., Ltd.	GFL-230
Microscope slide	Servicebio	
Upright electron microscope	JAPAN NIKON	NIKON ECLIPSE E100
Imaging system	JAPAN NIKON	NIKON DS-U3

2. Main experiment reagents

Reagent name	Manufacturer	Item No.
Xylene	Sinopharm Chemical Reagent Co., Ltd.	10023418
Gram staining solution kit	Servicebio	G1065
Neutral balsam	Sinopharm Chemical Reagent Co., Ltd.	10004160

II. Experimental procedure

1. Deparaffinize the paraffin sections to water: Place the sections in sequence in Xylene I for 20min - Xylene II for 20min - Anhydrous ethanol I for 5min - Anhydrous ethanol II for 5min - 75% Ethyl alcohol for 5min , wash with tap water.

2. Gram staining: Drop Gram staining solution A on the sections to stain for 10s-30s, wash and shake to make dry. Then add Gram staining solution B for 1min, wash and shake to make dry. And then take the Gram staining solution C with Pasteur pipette and rinse the tissue from the frosted surface side of the slide for differentiation, until the Gram staining solution C flow

down from the section is colorless, wash it and shake to get it dry. Dip the sections in Gram staining solution D for 1s and wash, microscopic examination, bake by oven to dry.

3. Dehydration: Quickly put the dried sections into three Jars of Anhydrous ethanol to dehydrate for 1s, 3s and 5s in sequence.

4. Transparency and sealing: Put the sections into clean Xylene for 5min for transparency, and then seal with neutral balsam.

5. Microscope inspection, image acquisition and analysis.

III. Interpretation of results

Gram-positive bacteria are purple, and Gram-negative bacteria are red.

IV. Precautions

1. The differentiation of Gram staining solution C in the staining process is a key step, and attention should be paid to the degree of differentiation.
2. Try not to repeat the staining. The differentiation solution removes the color of sections, and meanwhile also changes the permeability of the cell wall, which affects the accuracy of the final result.
3. When the temperature is low in winter, the staining time can be extended properly.

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Wuhan Servicebio Technology Co., Ltd.