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Golgi-cox Staining Experiment Report

I. Experimental Equipment and Reagents

1. Experimental equipment

Equipment name	Manufacturer	Model No.
Freezing microtome	Thermo	CRYOSTAR NX50
Slicing knife	Shanghai Leica Instrument Co., Ltd.	LEICA 819
Microscope slide	Citotest Labware Manufacturing Co., Ltd.	10127105P-G
Upright electron microscope	JAPAN NIKON	NIKON ECLIPSE E100
Imaging system	JAPAN NIKON	Nikon DS-U3
Digital slide scanner	3D HISTECH	Pannoramic 250

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Add: 5th Floor, C22 Building, Biopark, No. 388 Gaoxin 2nd Road, East Lake High-tech Developing Zone, Wuhan, Hubei, China, 430079 Website: <http://www.servicebio.com> Tel: 027-5111-3188

Gelatin potassium chromium sulfate mixture (It can be stored for 3 months under 4°C condition.)

Production of gelatin slides: Heat the gelatin potassium chromium sulfate mixture to liquid at 37°C (do not shake the solution, if small bubbles are formed in the solution, many small particles will remain on the slide during the subsequent production process), immerse ordinary slides into the liquid for a few seconds and take out, place on an iron rack and dry in a 37°C oven, store in a dry, cool and ventilated place.

1. Draw materials: After killing the rat, immediately take the brain tissue and fix it in the fixative for more than 48h.

2. Golgi-cox staining: Cut the rat brain tissue into tissue blocks with thickness of 2-3mm according to the tissue site to be observed, gently rinse the brain tissue with normal saline for several times, then place in a 45ml round bottom EP tube, add Golgi-cox staining solution to completely submerge the brain tissue, and place in a cool and ventilated place and avoid light to treat for 14 days (after 48h soak, change the new staining solution, and then change the new staining solution every 3 days, total of 14 days).

3. Tissue treatment: Immerse in distilled water for 3 times, pour in 80% glacial acetic acid to immerse the tissue overnight, wait for the tissue to become soft, then wash with distilled water, and place into 30% sucrose.

4. Sectioning: Cut the tissue into 100 microns with oscillating microtome, paste it on a gelatin slide, and dry in the dark overnight.

5. Developing and fixing: Treat dried tissue slides with concentrated ammonia water for 15 minutes, wash with distilled water for 1 minute, treat with acid hardening fixing solution for 15 minutes, wash with distilled water for 3 minute, dry and seal the section with glycerin gelatin.

6. Microscope inspection, image acquisition and analysis, obtain panoramic images of brain tissue by panoramic multi-layer scanning with digital slice scanner.

III. Interpretation of Results:

Neurons and dendritic spines are black, and the background is buff.

IV. Precautions:

1. The Golgi-cox staining solution is highly toxic, so the laboratory staff need to wear lab coat and disposable gloves to perform the experiment.
2. Each step during staining needs to be strictly protected from light.
3. The Golgi-cox staining solution should be stored in a cool, ventilated and dark place at room temperature. The fixing solution should be stored in a dark and ventilated place at 18°C-20°C.