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Glycine Silver Immersion Plating Nerve Staining

I. Experimental equipment and reagents

1. Experimental equipment

Equipment name	Manufacturer	Model No.
Dehydrator	DIAPATH	Donatello
Embedding center	Wuhan Junjie Electronics Co., Ltd.	JB-P5
Pathological microtome	Shanghai Leica Instrument Co., Ltd.	RM2016
Cooling plate	Wuhan Junjie Electronics Co., Ltd.	JB-L5
Tissue spreading water bath	Zhejiang Jinhua Kedi Instrumental Equipment Co., Ltd.	KD-P
Oven	Tianjin Leibo Terry Equipment Co., Ltd.	GFL-230
Microscope slide	Servicebio	
Upright electron microscope	JAPAN NIKON	NIKON ECLIPSE E100
Imaging system	JAPAN NIKON	NIKON DS-U3

2. Main experiment reagents

Reagent name	Manufacturer	Item No.
Anhydrous ethanol	Sinopharm Chemical Reagent Co., Ltd.	100092683
Xylene	Sinopharm Chemical Reagent Co., Ltd.	10023418
Glycine silver staining solution kit	Servicebio	G1052
Neutral balsam	Sinopharm Chemical Reagent Co., Ltd.	10004160

II. Experimental procedure

1. Deparaffinize the paraffin sections to water: Place the sections in sequence in Xylene I for 20min - Xylene II for 20min - Anhydrous ethanol I for 5min - Anhydrous ethanol II for 5min - 75% Ethyl alcohol for 5min , wash with tap water, wash with distilled water for several times, and wash clean.

2. Staining with Glycine silver staining solution C: Put the sections into the Glycine silver

staining solution C for 5min, wash with distilled water for 3 times.

3. Treatment with Glycine silver staining solution B: Put the sections into Glycine silver staining solution B for 3-5 minutes.

4. Staining with Glycine silver staining solution A: Take out the sections, quickly shake off the residual glycine silver staining solution B on the tissue, put it into glycine silver staining solution AI (preheated well at 45°C in advance), and observe the reduction effect at all times. Take out the sections after a few seconds and quickly put into the glycine silver staining solution A II (preheated at 45°C in advance). Quickly take out after a few seconds and wash with distilled water.

If the background of the sections is too dark, you can quickly treat it with glycine silver staining solution D for 1s, and wash with distilled water for 3 times.

5. Dehydration and sealing: Put the sections in order into Anhydrous ethanol I for 5min - Anhydrous ethanol II for 5min - Anhydrous ethanol III for 5min - Xylene I for 5min - Xylene II for 5min for transparency, at last, seal with neutral balsam.

6. Microscope inspection, image acquisition and analysis.

III. Interpretation of results:

Nerve axons, fibrils in nerve cells, neurofibrillary tangles, senile plaques and dendrites are black, background is golden yellow.

IV. Precautions:

1. It must be cleaned with distilled water before and after treatment with glycine silver solution, otherwise impurities will easily appear.
2. Only one section can be stained at a time, thus it's easy to control the effect, and the reduction process must be rapid, the film is easy to be dirty if it is not controlled well.
3. All staining solutions need to be protected from light. The glycine silver solution can be used 3-5 times in the same day, the reducing solution should be renewed each time.

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Wuhan Servicebio Technology Co., Ltd.