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Add: 5th Floor, C22 Building, Biopark, No. 388 Gaoxin 2nd Road, East Lake High-tech Developing Zone, Wuhan, Hubei, China, 430079 Website: <http://www.servicebio.com> Tel: 027-5111-3188

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Add: 5th Floor, C22 Building, Biopark, No. 388 Gaoxin 2nd Road, East Lake High-tech Developing Zone,
Wuhan, Hubei, China, 430079 Website: <http://www.servicebio.com> Tel: 027-5111-3188

Gross Oil Red Staining Experiment Report

I. Experimental equipment and reagents

1. Experimental equipment

Equipment name	Manufacturer	Model No.
Dissecting forceps		
Dissecting scissors		
Camera	Canon	D70

2. Main experimental reagents

Reagent name	Manufacturer	Item No.
Fixative	Servicebio	G1101
Oil red staining solution	Servicebio	G1015
Isopropyl alcohol	Sinopharm Chemical Reagent Co., Ltd.	80109218

II. Experimental procedure

1. Fixation: Take the blood vessels in target area, remove the adipose tissue around the blood vessels as much as possible with forceps, put the blood vessels in fixative to fix for more than 24h, take out the blood vessels tissues from the fixative, and wash twice with PBS.

2. Dissection: Remove the fat on the periphery of blood vessel when taking materials and carefully dissect the blood vessel longitudinally along the vessel wall with dissecting scissors.

3. Staining: Wash the well cut blood vessels with tap water for 5s slightly. Immerse the blood vessels in 60% isopropyl alcohol for 3s, and then immerse in the Oil red O staining solution. Stain in the staining solution at 37° in the dark for 60 minutes, then take out the blood vessels with forceps and immerse in 60% isopropyl alcohol for differentiation, the differentiation time began to differentiate at 1min until the fatty plaque in the vessels lumen is orange-red or bright red, and the other parts are nearly colorless, then wash with distilled water to terminate the differentiation.

4. Photographing: Take out the blood vessels, absorb excess water with filter paper, take a microscope slide and place it on a black or white background plate with scale, spread the blood vessel on the slide, select a place with good light, adjust the focus and exposure of the camera

well, then take photos. The scale on the background plate should also be taken into the photos.

III. Interpretation of results:

The lipid droplets are orange-red or bright red, and the other parts are nearly colorless.

IV. Precautions:

1. The staining results cannot be stored for a long time. It should be observed and take photos as soon as possible to save the experimental results.

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Wuhan Servicebio Technology Co., Ltd.