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Mucus Layer Preparation and PAS Staining

Experimental Report

I. Experimental Equipment and Reagents

1. Experimental equipment

Equipment name	Manufacturer	Model No.
Embedding center	Wuhan Junjie Electronics Co., Ltd.	JB-P5
Pathological microtome	Shanghai Leica Instrument Co., Ltd.	RM2016
Cooling plate	Wuhan Junjie Electronics Co., Ltd.	JB-L5
Tissue spreading water bath	Zhejiang Jinhua Kedi Instrumental Equipment Co., Ltd.	KD-P
Oven	Tianjin Leibo Terry Equipment Co., Ltd.	GFL-230
Microscope slide	Servicebio	G6012-1
Upright electron microscope	JAPAN NIKON	NIKON ECLIPSE E100
Imaging system	JAPAN NIKON	NIKON DS-U3

2. Main experiment reagents

Reagent name	Manufacturer	Item No.
Carnoy's Fluid	Servicebio	G1120
Methyl alcohol	Sinopharm Chemical Reagent Co., Ltd.	100141580
Anhydrous ethanol	Sinopharm Chemical Reagent Co., Ltd.	100092683
Xylene	Sinopharm Chemical Reagent Co., Ltd.	10023418
PAS Staining solution suit	Servicebio	G1008
Neutral balsam	Sinopharm Chemical Reagent Co., Ltd.	10004160

II. Experiment steps

1. **Draw materials and Fixation of intestinal tissue:** After the animal is killed, quickly cut out intestinal segment of target area (greater than 5 mm required). Keep the intestinal internal

contents when taking the intestinal segment samples. Do not rinse and squeeze the intestinal contents. Otherwise it will damage the mucus layer and affect the integrity of the tissue structure. After taking the sample, immediately put it into Carnoy's Fluid for fixation, for a minimum of 4 hours and no more than 7 days, thus to better preserve the intestinal mucus layer.

2. Tissue dehydration, embedding, sectioning: After the tissue is fixed, put into Methyl alcohol for 2 times, 40-60min each time, Anhydrous ethanol twice, 40-60min each time, Xylene transparency for 20-40min, Paraffin embedding, cut 4um paraffin sections, and bake the sections at 65 ° for an hour.

3. Deparaffinize the paraffin sections to water: Put the sections in order in Xylene I for 20min - Xylene II for 20min - Anhydrous ethanol I for 5min - Anhydrous ethanol II for 5min - 75% Alcohol for 5min, wash with tap water.

4. Put the sections into PAS staining solution B to stain for 10-15 minutes, wash with tap water and wash with distilled water twice.

5. Dip the sections into PAS staining solution A to stain for 25-30min and protect from light, rinse with running water for 5min.

6. Put the sections into PAS staining solution C for 30S, wash with tap water, differentiate with hydrochloric acid aqueous solution, wash with tap water, put in ammonia to return blue, and rinse with running water.

7. Dehydration and sections sealing : Put the sections in sequence into Anhydrous ethanol I for 5min - Anhydrous ethanol II for 5min - Anhydrous ethanol III for 5min - Xylene for 5min - Xylene II for 5min for transparency, seal the sections with neutral balsam.

8. Microscope inspection, image acquisition and analysis.

III. Interpretation of results:

The purple-red mucus layer is located between the intestinal contents and the intestinal mucosal epithelium, either is wrapped in the intestinal internal contents or is attached to the mucosal epithelium.

IV. Precautions:

1. Because the intestinal mucus is secreted by mucus cells into the intestinal cavity and attached to the mucosal epithelium, it is easy to fall off. And it is difficult to retain a very complete mucus layer. And the mucus belongs to carbohydrates, it must be fixed with water-free Carnoy's Fluid, otherwise the mucus layer will be dissolved.
2. PAS staining solution B should be replaced frequently, to avoid affecting the staining effect.
3. Avoid light when staining with PAS staining solution A.

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Wuhan Servicebio Technology Co., Ltd.