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Add: 5th Floor, C22 Building, Biopark, No. 388 Gaoxin 2nd Road, East Lake High-tech Developing Zone, Wuhan, Hubei, China, 430079 Website: <http://www.servicebio.com> Tel: 027-5111-3188

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Alcian Blue and Alizarin Red Staining Experimental Report

1. Experimental Equipment and Reagents

1.1 Experimental equipment

Name	Manufacturer	Model
Anatomic forceps		
Dissecting scissors		
Camera	Canon	EOS70D
Culture dish		
Specimen bottle	Servicebio	G6008-40ML

1.2 Main experimental reagents

Reagent name	Manufacturer	Article number
Alcian Blue	Sigma	A5268
Alizarin red	West Asia	
Glacial acetic acid	Sinopharm Group Chemical Reagent Co. Ltd	10000218
Ethanol	Sinopharm Group Chemical Reagent Co. Ltd	100092683
Glycerol	Sinopharm Group Chemical Reagent Co. Ltd	10010618
Potassium hydroxide	Sinopharm Group Chemical Reagent Co. Ltd	10017018
Acetone	Sinopharm Group Chemical Reagent Co. Ltd	10000428

2. Reagent Preparation

2.1 Alcian blue 8GS: 0.3 g alcian blue powder was diluted to 100 mL with 70% ethanol as solvent.

2.2 Alizarin red: 0.1g alizarin red powder was diluted to 100mL with 95% ethanol.

2.3 Alcian blue-alizarin red dye solution: mix 0.3% alcian blue: 0.1% alizarin red: glacial acetic acid: 70% ethanol in the ratio of 1:1:1:17.

2.4 1% glycerin:1% KOH in the ratio of 1:4

2.5 Preservation solution: 50% glycerin.

2.6 1% KOH: 1 g KOH powder is dissolved in 100 mL pure water.

3. Experimental Steps

3.1 Separation and fixation of skeletons: acetone treatment for 2~7 days (usually directly treated for 2 days during 13~16 days of gestation, and the skin and viscera of fetal rats were peeled off with two tweezers after pregnancy for more than 16 days), the tissues were completely immersed in acetone, and the treatment time could be appropriately prolonged according to the tissue size.

3.2 Acetone treatment for 2~7 days (usually direct treatment for 2 days at 13~16 days of pregnancy, and 7 days after peeling the skin and viscera of fetal rats with two tweezers), the tissue was completely immersed in acetone, and the treatment time can be appropriately prolonged according to the size of the tissue.

3.3 Take out the sample, put it in 20% glycerol and transparent for 1 week, until the skeleton and muscle are clear, and finally transfer to 50% glycerol to take pictures.

3.4 The sample was preserved in 50% glycerol.

4. The Results Were As Follows

4.1 The skeleton of mice is clear and transparent, the skeleton is purplish red, and the keratinized skin and muscle are blue, and the contrast is clear.

5. Precautions

5.1 Dehydration and acetone treatment time should not be too long, samples will appear uneven dyeing, especially claws and tail.

5.2 The dyeing time should be strictly controlled. If the dyeing time is too long, the whole skeleton will turn blue, if it is too light, the skeleton will fade easily.

5.3 1% KOH treatment should not be too long, otherwise the skeleton will disintegrate easily.

6. Pictures