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Add: 5th Floor, C22 Building, Biopark, No. 388 Gaoxin 2nd Road, East Lake High-tech Developing Zone, Wuhan, Hubei, China, 430079 Website: <http://www.servicebio.com> Tel: 027-5111-3188

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TTC Staining (Plant) Experimental Report

1. Experimental Equipment and Reagents

1.1 Experimental equipment

Name	Manufacturer	Model
Constant temperature water bath	Jiangyan Tianli Medical Instrument Factory	TL-42OD
Specimen bottle	Servicebio	G6008-40ML
Blade	Shanghai Leica Instrument Co., Ltd	LEICA 819
Dissecting forceps		
Camera	Canon	70D

1.2 Main experimental reagents

Reagent name	Manufacturer	Article number
TTC powder	Sigma	T8877
Fixative	Servicebio	G1101-500ML

2. Reagent Preparation

TTC solution: 1.0 g TTC powder dissolved in 1000 mL pure water, and mixed well, in reserve.

3. Experiment Steps

3.1 Soak Seeds in a specimen bottle containing 20mL tap water, and place them in a water bath at 35°C for 2-6h, or extend the soaking time until fully imbibition.

3.2 .The specimen bottle containing 20 mL TTC dye solution was heated in a water bath at 35°C in the dark for 30 min in advance. Bring out the seeds with tweezers, peel off the seed shell or seed coat, split the seed embryo into two parts, and place them in a TTC dyeing solution in a water bath at 37°C and stain for 0.5-1 h in the dark. After 30 min, observing the coloring of the tissue every 10 minutes, until the red of seed embryos is no longer darker or more.

3.3 After staining has been completed, clamp the tissue with tweezers, and empty the TTC solution in the sample bottle, and then pour 20mL fixative solution along the bottle wall. Finally, wrap the sample bottle in foil to store it in the dark at normal temperature for 24 hours, and then take photos to save the results.

4. The Results Were As Follows

The normal area of tissue was red, while the infracted area was pale.

5. Precautions

5.1 TTC dyeing solution needs to be stored at 4°C and away from light. It can be stored for one month, but should be recovered to normal temperature before usage.

5.2 It is necessary to avoid light during the dyeing process. The dyeing solution is easy to decompose when exposed to light, and the color of the dyed seeds will also fade when exposed to light.

5.3 Take pictures as soon as possible after staining to save experimental results.

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Wuhan Servicebio Technology Co., Ltd.