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Grimelins Silver Nitrate Staining Experiment Report

I. Experimental Equipment and Reagents

1. Experimental equipment

Equipment name	Manufacturer	Model No.
Dehydrator	DIAPATH	Donatello
Embedding center	Wuhan Junjie Electronics Co., Ltd.	JB-P5
Pathological microtome	Shanghai Leica Instrument Co., Ltd.	RM2016
Cooling plate	Wuhan Junjie Electronics Co., Ltd.	JB-L5
Tissue spreading water bath	Zhejiang Jinhua Kedi Instrumental Equipment Co., Ltd.	KD-P
Oven	Tianjin Leibo Terry Equipment Co., Ltd.	GFL-230
Microscope slide	Servicebio	
Upright electron microscope	JAPAN NIKON	NIKON ECLIPSE E100
Imaging system	JAPAN NIKON	NIKON DS-U3

2. Main experiment reagents

Reagent name	Manufacturer	Item No.
Anhydrous ethanol	Sinopharm Chemical Reagent Co., Ltd.	100092683
Xylene	Sinopharm Chemical Reagent Co., Ltd.	10023418
Silver nitrate	Sinopharm Chemical Reagent Co., Ltd.	10018461
Glacial acetic acid	Sinopharm Chemical Reagent Co., Ltd.	10000218
Sodium acetate	Sinopharm Chemical Reagent Co., Ltd.	10018718
Hydroquinone	Sinopharm Chemical Reagent Co., Ltd.	10011317
Anhydrous Sodium Sulfite	Sinopharm Chemical Reagent Co., Ltd.	S112302
Neutral balsam	Sinopharm Chemical Reagent Co., Ltd.	10004160

II. Dyeing Solution Preparation

1% silver nitrate solution: 1g Silver nitrate + 100mL ultrapure water

Acetate buffer (mother liquid) : Solution A: 1.21ml glacial acetic acid + 98.7ml ultrapure water

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Add: 5th Floor, C22 Building, Biopark, No. 388 Gaoxin 2nd Road, East Lake High-tech Developing Zone, Wuhan, Hubei, China, 430079 Website: <http://www.servicebio.com> Tel: 027-5111-3188

Solution B: 2.72g sodium acetate + 100ml ultrapure water

Reducing solution: 1g Hydroquinone + 2.5g anhydrous sodium sulfite + 100ml ultra-pure water

III. Experiment Steps

- 1. Deparaffinize the paraffin sections to water:** Place the sections in order in Xylene I for 20min - Xylene II for 20min - Anhydrous ethanol I for 5min - Anhydrous ethanol II for 5min - 75% Ethyl alcohol for 5min, wash with tap water, and 3-5 times with distilled water.
- 2. Prepare silver solution incubation medium:** Mix solution A: solution B=1:9 into acetate buffer, take 3ml 1% silver nitrate + 10ml acetate buffer + 87ml ultrapure water.
- 3. Silver nitrate staining:** Place the silver solution incubation medium at 60 °C for preheating, and immerse the sections in the silver solution for incubation at 60 °C for 3 hours.
- 4. Reducing dyeing:** Take out the sections from the silver solution, quickly dip it in distilled water for 1 time, then put the sections into the 45°C preheated reducing solution for 1 hour, take out the sections, dip and wash in distilled water for several times.
- 5. Dehydration and transparency:** Put the sections in order into Anhydrous ethanol I for 5min - Anhydrous ethanol II for 5min - Anhydrous ethanol III for 5min - Xylene I for 5min - Xylene II for 5min for transparency, at last, seal with neutral balsam.
- 6. Microscope inspection,** image acquisition and analysis.

IV. Interpretation of Results:

The argyrophil cells are mainly in the shape of vases, spindle, tadpole, circle and badminton racket. The nucleus is non-staining, the intracytoplasmic argyrophilic grain are plentiful and in brown-black color.

V. Precautions:

1. Both the mixed solution and the silver solution incubation medium should be mixed freshly and use immediately, the mother liquid of the mixed solution can be stored for a long time.
2. During the dyeing process, rapid dip and wash in distilled water requires quick operation.