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1.2 Major reagents

Name	Producer	Code
Ethanol	SCRC	100092683
Xylene	SCRC	10023418
EVG dye solution set	Servicebio	G1042
Neutral gum	SCRC	10004160

2 Procedure

2.1 Dewaxing as followed:

Xylene I for 20 min;
Xylene II for 20 min;
100% ethanol I for 5 min;
100% ethanol II for 5 min;
75% ethanol for 5 min;
Rinsing with tap water ;

2.2 EVG A, EVG B and EVG C were prepared into EVG solution according 2 days in advance to the ratio of 5:2:2. Then stain with EVG B for 5 min, rinsing with tap water.

2.3 Treat the section repeatedly with double diluted EVG B and tap water until the background was separated. Control the degree of differentiation under a microscope, until the elastic fibers are purple-black, and the background is gray-white and almost colorless.

- 2.4 VG dye solution was prepared according to the ratio of 9ml EVG E and 1ml EVG D. Configure according to actual usage. Stain with VG solution for 1-3 min. Then washing it fast with water, dehydration quickly with three cup of anhydrous ethanol. It should be noted that the staining time depends on the composition of the elastic fibers in the tissue. If the staining time is too short, the collagen will be light in color, and if the staining time is too long, the elastic fibers will fade.
- 2.5 Xylene I for 20 s; Xylene II for 5 min; Finally seal with neutral gum. Please do not use xylene used in other experiments.
- 2.6 Observe with microscope inspection, image acquisition and analysis.

3 Results

Color	Result
Purple-black	Elastic fibers
Red	Collagen fibers
Yellow	Background

4 Precautions

- 4.1 During the differentiation, the elastic fibers should be in the form of purple-black filaments. Do not over-differentiate, otherwise the elastic fibers will fade. If the differentiation is insufficient, the effect of VG counterstaining will be poor.

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Wuhan Servicebio Technology Co., Ltd.