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DAB Staining Experimental Report

1. Experimental Equipment and Reagents

1.1 Experimental equipment

Name	Manufacturer	Model
Clean workbench	AIRTECH	SW-CJ-1FD
Specimen bottle	Servicebio	
Anatomy tweezers		
Camera	Canon	70D

1.2 Main experimental reagents

Reagent name	Manufacturer	Article number
DAB dye solution (plant)	Servicebio	G1022
Absolute ethanol	Sinopharm Group Chemical Reagent Co. LTD	100092683

2. Experimental Steps

2.1 Collect seedlings or root tips of plants, and wash them with pure water. Then put them on filter paper with tweezers to dry out, and immerse them into a 10mL round-bottom EP tube containing 5mL DAB dye solution A, dipping for 2-6h at normal temperature in the dark.

2.2 Take out the plant seedlings or leaves with tweezers, and immerse them in pure water and rinse 5 times, Then put them on the filter paper to dry out, and soak them into 95% ethanol at 40°C for 3-16h.

2.3 Take out the plant seedlings or leaves with tweezers, and soak them in pure water again and rinse 5 times. Then put them on the filter paper to dry out and take pictures. If you need to save the sample, store the stained sample in a sample bottle containing 20 mL DAB dye solution B at room temperature.

3. The Results Were As Follows

The parts that produce hydrogen peroxide in the root tips or leaves of the plant would be dark brown, while the others were light brown or colorless, or in the color of the plant pigment itself.