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Genomic DNA Extraction Report

1 Laboratory Equipment and Reagents

1.1 Laboratory equipment

Equipment	Manufacturers	Model
Centrifuge	DragonLab	D3024R
Clean bench	Suzhou Antai Air Tech Co., Ltd	SW-CJ-1FD
Standard reagent type pure water meter	Qingdao Fulum Technology Co., Ltd	FBZ2001-up-p
Ultramicro spectrophotometer	Thermo	NanoDrop2000

1.2 Reagents

Reagents	Manufacturers	Order
(Tissue/Cell/Blood) Genomic Extractin DNA Kit	Servicebio	G3633
Isopropyl alcohol	Sinopharm Group Chemical Reagent Co., Ltd.	80109218
Centrifuge tube	Servicebio	
TIP	Servicebio	

2 Genomic DNA Extraction Experimental steps

2.1 Samples preparation

2.1.1 For blood, please use 50-100 μ L fresh, frozen or anticoagulant adding blood. If less than 200 μ L, please make up with **buffer GA** to 200 μ L.

2.1.2 If the sample is blood from poultry, birds, amphibians, of which red blood cells have nucleolus, the amount should be reduced to 5-20 μ L and adjust the volume to 200 μ L with **buffer GA**.

2.1.3 The adherent cells should be treated to cell suspension first, then centrifuge the cells for 5 min at 1,000 rpm, then discard the flow-through and re-suspend cell pellet in 200 μ L **buffer GA**.

2.1.4 Animal tissue (2-30mg) should be treated to cell suspension first, then centrifuge the cells for 1 min at 12,000 rpm, then discard the flow-through and re-suspend cell pellet in 200 μ L **buffer GA**.

2.2 Add 20 μ L **Proteinase K Stock Solution**, mix thoroughly by vortex. If the sample is tissue: incubate at 56°C for 30-60 min until the tissue is completely lysed (mix thoroughly by vortex per 10 min).

2.3 Add 200 μ L **Buffer GB** to the sample, mix thoroughly for 6-8 times, and incubate at 70°C for

10 min to yield a homogeneous solution. Briefly centrifuge the 1.5 ml microcentrifuge tube to remove drops from the inside of the lid.

2.4 Add 200 μ L ethanol to the sample, and mix thoroughly by vortex for 15 s. A white precipitate may form on addition of ethanol. Briefly centrifuge the 1.5 ml microcentrifuge tube to remove drops from the inside of the lid.

2.5 Pipet the mixture from step 4 into the Bind DNA Mini Columns (in a 2 ml collection tube) and centrifuge at 10000g for 30 s. Discard flow-through and place the spin column into the collection tube.

2.6 Add **500 μ L Buffer PD** to Bind DNA Mini Columns, and centrifuge at 10000g for 1 min , then discard the flow-through and place the spin column into the collection tube.

2.7 Add **600 μ L Buffer PW** (Ensure ethanol has been added) to Bind DNA Mini Columns, and centrifuge at 10000g for 1 min. Discard the flow-through and place the spin column into the collection tube.

2.8 Repeat Step 7.

2.9 Centrifuge at 10000g for 2 min to dry the membrane completely. Note: The residual ethanol of buffer PW may have some effect in downstream application.

2.10. Place the Bind DNA Mini Column in a new clean 1.5 ml microcentrifuge tube, and pipet 50-100 μ L **Elution Buffer or ddH₂O** directly to the center of the membrane. Incubate at room temperature (15-25°C) for 2-5 min, and then centrifuge for 2 min at 10000g.