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### 2、Reagents

| Reagents                                          | Company            | Catlog      |
|---------------------------------------------------|--------------------|-------------|
| RIPA buffer                                       | Servicebio         | G2002       |
| 50* cocktail                                      | Servicebio         | G2006       |
| PMSF (100mM)                                      | Servicebio         | G2008       |
| Phosphatase inhibitor                             | Servicebio         | G2007       |
| BCA protein quantitative detection kit            | Servicebio         | G2026       |
| 5*loading buffer                                  | Servicebio         | G2013       |
| SDS-PAGE kits                                     | Servicebio         | G2003       |
| Protein Marker                                    | Thermo (Fermentas) | 26617       |
| Electrophoretic buffer                            | Servicebio         | G2018       |
| Coomassie brilliant blue dye                      | Servicebio         | G2021-250ML |
| Decolorizing solution of Coomassie Brilliant blue | Servicebio         | G2022-250ML |

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Add: 5th Floor, C22 Building, Biopark, No. 388 Gaoxin 2nd Road, East Lake High-tech Developing Zone, Wuhan, Hubei, China, 430079 Website: <http://www.servicebio.com> Tel: 027-5111-3188

### (1) Sample Preparation (Total protein)

**1.1 For suspension cells:** Centrifuge at 2000rpm, 4°C for 5 min to collect cells, add RIPA buffer. (250μL RIPA buffer/10<sup>6</sup> cells (Protease and phosphatase inhibitors should be added before use) . If you wish to increase protein concentration, then add less lysis buffer.

#### 1.2 For adherent cells

1.2.1 Place the cell culture dish on ice and wash the cells three times with ice-cold PBS.

1.2.2 Aspirate the PBS, then add ice-cold RIPA buffer (Protease and phosphatase inhibitors should be added before use).

1.2.3 Scrape adherent cells off the dish , then gently transfer the cell suspension into a 1.5mL microcentrifuge tube.

1.3 Maintain constant agitation for 30 min on ice bath, use pipette to mix it up.

1.4 Centrifuge at 12,000 rpm , 4°C for 10 min. Gently remove the tubes from the centrifuge and place on ice, transfer the supernatant to a fresh tube.

### (2) Preparation of lysate from tissues

2.1 Wash the sample with cold PBS for twice or three times , then homogenize with an electric homogenizer. Volumes of lysis buffer must be determined in relation to the amount of tissue.

2.2 Maintain constant agitation for 30 min on ice bath, use pipette to mix it up.

2.3 Centrifuge at 12,000 rpm ,4°C for 10 min . Gently aspirate the supernatant and place in a fresh tube kept on ice.

### (3) SDS-PAGE

3.1 Washing glass plates.

3.2 Preparing and loading the gel.

3.3 Align the glass plate and put it into the gel maker to align and clamp tightly to avoid gel leakage

3.4 Different concentrations of separation gel were prepared according to the experimental requirements, and mixed evenly immediately after adding TEMED, and the separation gel was poured to an appropriate height. After about 45min, pour the water from the upper layer of the separation glue and drain the remaining water with absorbent paper.

| Reagents                     | 8%   | 10% | 12% | 15% | 18% | 20%  |
|------------------------------|------|-----|-----|-----|-----|------|
| H <sub>2</sub> O mL          | 4.63 | 4   | 3.3 | 2.3 | 1.3 | 0.63 |
| 30% acrylamide (29: 1)<br>mL | 2.67 | 3.3 | 4   | 5   | 6   | 6.67 |
| 1.5M TRIS—HCl(PH 8.8)<br>mL  | 2.5  | 2.5 | 2.5 | 2.5 | 2.5 | 2.5  |

|               |       |     |     |     |     |     |
|---------------|-------|-----|-----|-----|-----|-----|
| 10%SDS mL     | 0.1   | 0.1 | 0.1 | 0.1 | 0.1 | 0.1 |
| AP mL         | 0.1   | 0.1 | 0.1 | 0.1 | 0.1 | 0.1 |
| TEMED $\mu$ L | 5     | 5   | 5   | 5   | 5   |     |
| Total         | 10 mL |     |     |     |     |     |

| Reagents                  | Concentration 5% |
|---------------------------|------------------|
| H <sub>2</sub> O mL       | 2                |
| 30% acrylamide (29: 1) mL | 0.5              |
| 1M TRIS-HCl(PH 6.8)mL     | 0.5              |
| 10%SDS $\mu$ L            | 40               |
| AP $\mu$ L                | 30               |
| TEMED $\mu$ L             | 4                |
| Total mL                  | 3                |

3.5 Prepare 5% concentrated gel according to the above formula. After adding TEMED, mix well immediately, and then the gel can be filled. Fill the remaining space with the concentrate and insert the comb into the concentrate. Make sure there are no bubbles under the comb.

3.6 After stacking gel has polymerized, remove comb carefully and begin the electrophoresis.

3.7 Add enough electrophoresis solution for sample electrophoresis. Concentrated gel voltage is 75V, and separation gel uses 120V. The electrophoresis was stopped when the bottom of bromophenol blue was about 1cm.

#### (4) Coomassie brilliant blue staining and subsequent experimental steps

4.1 Put the gel into a glass plate, add Coomassie brilliant staining solution, and incubate at 70 °C for 5 min, and then place it on a shaker for 90 min;

4.2 After recovery, wash the gel with water, add the decolorizing solution, and put it on the decolorizing shaker overnight;

4.3 Remove the decolorizing solution, clean the gel with water and take photos of the target band.

#### (5) Results and analysis (see excel gray analysis table for details)