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Add: 5th Floor, C22 Building, Biopark, No. 388 Gaoxin 2nd Road, East Lake High-tech Developing Zone, Wuhan, Hubei, China, 430079 Website: <http://www.servicebio.com> Tel: 027-5111-3188

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Luxol Fast Blue (LFB) Staining Experimental Report

1. Lab equipment and reagents

1.1 Lab equipment

Items	Manufacturer	Model
Dehydrator	DIAPATH	Donatello
embedding machine	Wuhan Junjie Electronics Co., Ltd.	JB-P5
Pathology microtome	Shanghai Leica Instruments Co., Ltd.	RM2016
Frozen platform	Wuhan Junjie Electronics Co., Ltd.	JB-L5
Water Bath-Slide Drier	Zhejiang Jinhua Kedi Instrumental Equipment CO.,LTD	KD-P
Laboratory oven	Tianjin Labotery Instrument Equipment Co., Ltd.	GFL-230
Microscope slide	Servicebio	
Upright optical microscope	Nikon Japan	Nikon Eclipse E100
Imaging system	Nikon Japan	NIKON DS-U3

1.2 Chemical Reagents

Items	Manufacturer	Model
Absolute alcohol	Sinopharm Chemical Reagent Co., Ltd.	100092683
Xylene	Sinopharm Chemical Reagent Co., Ltd.	10023418
Luxol fast blue staining kit	Servicebio	G1030
Hematoxylin-Eosin solution	Servicebio	G1002
Neutral balsam	Sinopharm Chemical Reagent Co., Ltd.	10004160

2. Experimental steps

2.1 Paraffin section deparaffinization and rehydration: put the slides into xylene I 20 minutes-xylene II 20 minutes-absolute ethanol I 5 min-absolute ethanol II 5 min-75% alcohol for 5 min, tap water washing.

2.2 Luxol fast bluestaining: placed luxol fast blue staining A is in a 60°C oven and preheated for 30 minutes. Put he slides intoluxol fast blue staining A and cover with membrane for 1 hour, then take out the slides quickly wash with tap water.

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Wuhan Servicebio Technology Co., Ltd.

2.3 Background differentiation: put the slides into luxol fast blue staining B for 2s and directly immerse in luxol fast blue staining C for differentiation 15s, washing to terminate differentiation, microscope examination, repeated differentiation, washing and microscopic examination, until the myelin sheath is blue, the background is almost colorless;

2.4 Eosin counterstaining: put the slides in the oven at 65 degrees to dry it (about 30mins until the slides are dry) and take out and cool it, put into 95% ethanol and Eosin counterstaining.

2.5 Dehydration and sealing: put the slides into absolute ethanol I 5min-anhydrous ethanol II 5min-anhydrous ethanol III 5min-xylene I 5min-xylene II 5min -transparent, neutral balsam sealing.

2.6 Microscope examination, images collection and analysis.

3. Results

Nerve myelin sheath is blue, and other components are almost colorless.

4. Note

1. After stained, do not wash them in water for a long time due to the color will be fade.
2. Pay attention to the degree of alternating differentiation, if over-differentiation, the color of the myelin sheath will be fade, insufficient differentiation, the boundaries are unclear.
3. If over-differentiation, you can re-stain the slides in the staining solution.