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Experimental Report on Waxy Oil Red Staining of Leaves

1. Experimental equipments and reagents

1.1 Experimental equipments

Name	Manufactor	Model
Frozen slicer	Thermo	CRYOSTAR NX50
Upright optical microscope	Nikon Japan	Nikon Eclipse E100
Imaging system	Nikon Japan	Nikon DS-U3

1.2 Main experimental reagents

Reagent name	Manufactor	Article number
sucrose	Sinopharm	
Oil red o saturated liquid	Servicebio	G1015
Glycerin gelatin sealed tablets	Servicebio	G1402
Isopropanol	Sinopharm	

2. Experimental steps

- 2.1 Preparation: The leaves are fixed by FAA and then dehydrated by sucrose, or fresh plant leaves, prepared frozen slices (longitudinal cut), and dried at room temperature for about 15min;
- 2.2 Staining: Stain into oil-saturated O liquid for 8-10 minutes (protect from light and cover);
- 2.3 Differentiation: Differentiate with two cylinders of 60% isopropyl alcohol for about 10 s;
- 2.4 Washing: dip in distilled water, move gently to prevent the displacement of fat droplets;
- 2.5 Sealing: filter paper to dry excess water around the tissue and seal with glycerin gelatin;
- 2.6 Photography: avoid light, store slices at room temperature, observe and take pictures with a microscope, 40x objective lens and oil lens are recommended for observation.

3. The results were as follows

The wax was orange-red.