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Leaf Fungal Hyphae WGA-PI Staining Experimental Report

1. Experimental equipment and reagents

1.1 Experimental equipment

Name	Manufactor	Model
CO ₂ incubator	SANYO	MCO-15AC
Clean workbench	Su Jing An Tai	SW-CJ-1FD
Slides	Servicebio	
Cover glass	Jiangsu Shitai	10127105P-G
Orthostatic microscope	Nikon, JAPAN	Nikon Eclipse E100
Image system	NIKON, JAPAN	Nikon DS-U3

1.2 Main experimental reagents

Reagent name	Manufactor	Article number
Carnot fixative	Servicebio	G1120
Potassium hydroxide	Sinopharm Group Chemical Reagent Co. LTD	
WGA-AF488 storage solution	Servicebio	
PBS	Servicebio	G0002
PI stock solution	Servicebio	
Anti-fluorescence quenching sealed tablets	Servicebio	G1401

2. Experimental steps

2.1 Dye preparation: Add 10 μ L WGA storage solution and 20 μ L PI stock solution to 970 μ L 0.2% Tween-PBS, and mix thoroughly.

2.2 Fix: Use Carnot fixative to fix sample (Remove chlorophyll at the same time).

2.3 Transfer the leaves into 10% KOH solution, and seal the tube with a parafilm to prevent the tube from collapsing, and incubate at 85°C for a certain period of time (Generally speaking, corn leaves are treated for 4h, and the barley leaves are treated for 2h.).

2.4 Staining pretreatment: Tissues are washed 4 or 5 times with PBS.

2.5 Staining: Add dye solution to the centrifuge tube containing blades, and filter 4 times with vacuum filter where each time lasts for 5min, and each interval is 5min under normal pressure.

2.6 Photographing: The leaves were washed 2 or 3 times with PBS, and sealed with anti-fluorescence quenching to be stored at 4°C in the dark. And photographed with a fluorescence microscope or a confocal microscope.

3. The results were as follows

Fungal hyphae shows green fluorescence, cell walls of plant show red fluorescence.