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Plant Softening Paraffin Embedding Experiment Report

1. Experimental equipment and reagents

1.1 Experimental equipment

Name	Manufacturer	Model
Vacuum tissue dehydrator	Thermo Fisher	Excelsior AS
Embedding machine	Wuhan Junjie Electronics Co., Ltd	JB-P5
Frozen platform	Wuhan Junjie Electronics Co., Ltd	JB-L5
Electric blast drying oven	Tianjin Bo Ruite Instrument Co., Ltd	GFL-230
Electronic scale	Electronic	SF-400

1.2 Main experimental reagents

Reagent	Manufacturer	Article number
Xylene	Sinopharm Group Chemical Reagent Co. Ltd	10023418
Ethanol	Sinopharm Group Chemical Reagent Co. Ltd	100092683
Ethylenediamine anhydrous	Tianjin Damao Chemical Reagent Factory	107-15-3
Plant softener	Wuhan Seville Biotechnology Co., Ltd	G1115

2. Plant softening process

2.1 According to the location and requirements of the material, it is best to thin the tissue to 2~3mm, and then put the printed with the tissue in the dehydration box

2.2 Put the dehydration box into the plant softening solution, seal it, and put it into a 55 °C blast drying box to soften. The replacement period of softening solution is 7 days, and the softening degree is observed once a week.

2.3 Rinse the softened plants with running water for 30 minutes, soak them in 15% ethanol for 2 hours, then transfer to a vacuum tissue dehydrator to dehydrate 30% alcohol for 1 hour, 50% alcohol for 1 hour, 75% alcohol for 1 hour, 85% alcohol for 1.5 hours, 90% alcohol for 1.5 hours, 95% alcohol for 1.5 hours, 100% alcohol for 1 hour, 100% alcohol for 1 hour, xylene 8 min, xylene 20 min, xylene 20 min, paraffin 30 min, paraffin for 1 hour, paraffin for 2 hours.

2.4 Put the melted wax into the embedding frame, and before the wax solidifies, remove the tissue from the dewatering box and put it into the embedding frame according to the requirements of the embedding surface and affix the corresponding label. The wax is cooled at -20 °C, and the wax blocks is taken after the wax is solidified

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Wuhan Servicebio Technology Co., Ltd.

2.5 The wax blocks were placed on a paraffin slicer with a thickness of 5 μ m. The slice floats on the 40 °C warm water of the spreading machine, pick up the slides and bake in the oven at 60 °C. And then take out and store at room temperature and set aside.