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Experimental Report on Plant Tissue Safranin o-Fast Greenstaining

1. Experimental equipment and reagents

1.1 Experimental equipment

Name	Manufacturers	Model
Dehydrator	DIAPATH	Donatello
Encapsulation machine	Wuhan Junjie Electronics Co., Ltd.	JB-P5
Pathological slicer	Shanghai Leica Instruments Co., Ltd.	RM2016
Frozen	Wuhan Junjie Electronics Co., Ltd.	JB-L5
Organizing a stall machine	Cody Equipment Co., Ltd., Jinhua City, Zhejiang Province	KD-P
Oven	Tianjin Leiboi Instruments and Equipment Co., Ltd.	GFL-230
Slides	Servicebio	
Positive optical microscope	Nikon, Japan	Nikon Eclipse E100
Imaging systems	Nikon, Japan	Nikon DS-U3

1.2 Main experimental reagents

Reagent Name	Manufactor	Article Number
Waterless ethanol	China Pharmaceutical Group Chemical Reagents Co., Ltd.	100092683
Xylene	China Pharmaceutical Group Chemical Reagents Co., Ltd.	10023418
BioDewax and Clear Solution	Servicebio	G1128
Safranin O-Fast Green staining solution	Servicebio	G1031
Neutral gum	China Pharmaceutical Group Chemical Reagents Co., Ltd.	10004160

2. Experimental steps

2.1 The sections were rehydrated in two changes of BioDewax and Clear Solution and 100%I% - 100%II - 75% alcohol. Each step takes anywhere for 5 minutes. Finally, rinse with running water.

2.2 Safranin O staining: Put sections into safranin O staining solution for 15-30 s and three cylinders of anhydrous ethanol for rapid dehydration.

2.3 Decolorization: put slides into 50%, 70% and 80% alcohol for 3~8 s.

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Add: 5th Floor, C22 Building, Biopark, No. 388 Gaoxin 2nd Road, East Lake High-tech Developing Zone,
Wuhan, Hubei, China, 430079 Website: <http://www.servicebio.com> Tel: 027-5111-3188

2.4 Fixed green staining: sections into plant solid green staining solution staining 6~20 s, anhydrous ethanol three-cylinder dehydration.

2.5 Put sections into three cylinders of xylene for 5 min. Finally, mount the tissue section with neutral balsam.

2.6 Observed under microscope and took images.

3. The results were as follows

The cytoderm which has lignified appear red and the color of cellulose cell wall is green.

4. Precautions

4.1 Wash quickly after stain safranin O, because it is easy to fade in water.

4.2 The film after dyeing should have a clear morphology and structure, and no excess solid green dyeing solution remains.