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Hard Plant Safranin o-Fast Green (Hard Plant) Stain Report

1. Experimental equipment and reagents

1.1 Equipment

Brand name	Manufactures	Model
Oven	Tianjin Labotery instrument equipment Co. Ltd	GFL-230
Microslide	Servicebio	
Optical Microscopy	NIKON	Nikon Eclipse E100
Imaging system	NIKON	Nikon DS-U3

1.2 Reagent

Reagent	Manufactures	Cat
Ethanol absolute	Sinopharm Chemical Reagent Co.,Ltd	100092683
Xylenes	Sinopharm Chemical Reagent Co.,Ltd	100023418
Safranin O-Fast Green staining solution(plant)	Servicebio	G1031
Neutral balsam	Sinopharm Chemical Reagent Co.,Ltd	10004160
Ethylene glycol ethyl ether acetate	Macklin	E808814-2.5L

2. Experiment procedure

2.1 Remove the embedding agent and rehydration

Put sections into ethylene glycol ethyl ether acetate I 6 h at 37 °C and ethylene ether acetate II overnight at 37 °C. Then put them into ethylene glycol ethyl ether acetate III 10-15 min minutes at room temperature and ethylene glycol ethyl ether acetateIV 10-15 min minutes at room temperature. The sections were rehydrated in 100%I% - 100%II - 95% - 90% - 80% alcohol. Each step takes anywhere for 10 minutes. Finally, rinse with running water.

2.2 Safranin O staining

Put sections into safranin O staining solution for 15-30 s and three cylinders of anhydrous ethanol for rapid dehydration.

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Add: 5th Floor, C22 Building, Biopark, No. 388 Gaoxin 2nd Road, East Lake High-tech Developing Zone, Wuhan, Hubei, China, 430079 Website: <http://www.servicebio.com> Tel: 027-5111-3188

Put sections into three cylinders of 50%, 70%, 80% alcohol respectively for 3-8 s. The concentrations of alcohol are 50,70 and 80 percent.

2.4 Fast Green staining

Put sections into fast green staining solution for 4-6 s and three cylinders of anhydrous ethanol for 5 s in turn.

2.5 Transparent and mount

Put sections into three cylinders of xylene for 5 min. Finally, mount the tissue section with neutral balsam.

2.6 Microscope observation and collection and analysis images.

3. Result determination

The cytoderm which has lignified appear red and the color of cellulose cell wall is green.

4. Matter need attention

1. Wash quickly after stain safranin O, because it is easy to fade in water.