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Tunel & Immunofluorescence Double Staining Protocol (Cell Climbing Slides)

1. Apparatus and reagents

1.1 Apparatus

Name	Producer	Model
Glass microscope slide	Servicebio	
Vortex	Servicebio	MX-F
Micro-centrifuge	Servicebio	D1008E
Rocker	Servicebio	TSY-B
Pipettor	Dragon	KE0003087/KA0056573
Liquid blocker pen	Servicebio	G6100
Ortho-Fluorescent Microscopy	Nikon	NIKON ECLIPSE C1
Imaging system	Nikon	NIKON DS-U3

1.2 Major reagents

Name	Producer	Model
BSA	Servicebio	G5001
PBS solution	Servicebio	G0002
Permeabilize solution	Servicebio	G1204
Tunel assay kit	Servicebio	G1501
Primary antibody		
Fluorescent-labelled secondary antibody		
DAPI	Servicebio	G1012
anti-fade mounting medium	Servicebio	G1401

2. Procedure

2.1 Permeabilization: dry the cell climbing slides slightly. Mark the objective area with liquid blocker pen, where add 50-100 μ l of permeabilize working solution. Incubate for 20 min at room temperature. Wash three times with PBS solution, 5 min each.

2.2 Equilibrium at room temperature: After the slices are slightly dried, buffer is added to the tissues in the circle, and the buffer is incubated at room temperature for 10 minutes

2.3 Tunnel reaction: Take appropriate amount of TDT enzyme, dUTP and buffer in the tunnel kit according to the number of slices and tissue size and mix at 1:5:50 ratio. Prepare this reaction solution according to demand before use. Add this mixture to objective tissue placed in a flat wet box, incubate at 37°C for 2 h. be sure to keep the wet box moist by adding water.

2.4 BSA blocking: wash slides three times with PBS (pH 7.4) in a Rocker device, 5 min each. Add 3% BSA to cover the marked tissue to block non-specific binding for 30 min BSA.

2.5 Primary antibody: throw away the blocking solution slightly. Incubate slides with primary antibody (diluted with PBS appropriately) overnight at 4 °C, placed in a wet box containing a little water.

2.6 Secondary antibody: wash slides three times with PBS (pH 7.4) in a Rocker device, 5 min each. Then throw away liquid slightly. Cover objective tissue with secondary antibody (appropriately respond to primary antibody in species), incubate at room temperature for 50 min in dark condition.

2.7 DAPI counterstain in nucleus: wash three times with PBS (pH 7.4) in a Rocker device, 5 min each. Then incubate with DAPI solution at room temperature for 10 min, kept in dark place.

2.8 Mount: wash three times with PBS (pH 7.4) in a Rocker device, 5 min each. Throw away liquid slightly, then coverslip with anti-fade mounting medium.

2.9 Microscopy detection and collect images by Fluorescent Microscopy. DAPI glows blue by UV excitation wavelength 330-380 nm and emission wavelength 420 nm; FITC glows green by excitation wavelength 465-495 nm and emission wavelength 515-555 nm; CY3 glows red by excitation wavelength 510-560 nm and emission wavelength 590 nm.

3. Results

The nuclei stained by DAPI are blue under the excitation of ultraviolet light, the tunnel kit is labeled with FITC fluorescein, the positive apoptotic cells are labeled with green light, and the primary antibody is labeled with red.