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(Red) Tunel Assay Protocol (Cell Climbing Slides)

1. Apparatus and reagents

1.1 Apparatus

Name	Producer	Model
Glass microscope slide	Servicebio	
Vortex	Servicebio	MX-F
Micro-centrifuge	Servicebio	D1008E
Rocker	Servicebio	TSY-B
Pipettor	Dragon	KE0003087/KA0056573
Liquid blocker pen	Servicebio	G6100
Ortho-Fluorescent Microscopy	Nikon	NIKON ECLIPSE C1
Imaging system	Nikon	NIKON DS-U3

1.2 Major reagents

Name	Producer	Model
PBS solution	Servicebio	G0002
Permeabilize solution	Servicebio	G1204
DAPI	Servicebio	G1012
anti-fade mounting medium	Servicebio	G1401
Tunel assay kit	Servicebio	G1502

2. Procedure

2.1 Permeabilization: dry the cell climbing slides slightly. Mark the objective area with liquid blocker pen, where add 50-100 μ l of permeabilize working solution. Incubate for 20 min at room temperature. Wash three times with PBS solution, 5 min each.

2.2 Equilibrium at room temperature: After the slices are slightly dried, buffer is added to the tissues in the circle, and the buffer is incubated at room temperature for 10 minutes

2.3 Tunel reaction: Take appropriate amount of TDT enzyme, dUTP and buffer in the tunel kit according to the number of slices and tissue size and mix at 1:5:50 ratio. Prepare this reaction solution according to demand before use. Add this mixture to objective tissue placed in a flat wet box, incubate at 37°C for 2 h. be sure to keep the wet box moist by adding water

2.4 DAPI counterstain in nucleus: wash three times with PBS (pH 7.4) in a Rocker device, 5 min each. Then incubate with DAPI solution at room temperature for 10 min, kept in dark place.

2.5 Mount: wash three times with PBS (pH 7.4) in a Rocker device, 5 min each. Throw away liquid slightly, then coverslip with anti-fade mounting medium.

2.6 Photomicrograph:

The sections were observed under a fluorescence microscope and images were collected. (DAPI ultraviolet excitation wavelength 330-380nm, emission wavelength 420nm, blue light emission; CY3 excitation wavelength 510-561nm, emission wavelength 590 nm, red light emission)

3. Results

The nuclei stained by DAPI are blue under the excitation of ultraviolet light, the tune3 kit is labeled with CY3 fluorescein, and the positive apoptotic nuclei are red

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Wuhan Servicebio Technology Co., Ltd.