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1.2 Major reagents

Name	Producer	Code
Ethanol	SCRC	100092183
Biodewax and Clear Solution	SERVICEBIO	G1128
PBS solution	Servicebio	G0002
Permeabilize solution	Servicebio	G1204

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Add: 5th Floor, C22 Building, Biopark, No. 388 Gaoxin 2nd Road, East Lake High-tech Developing Zone, Wuhan, Hubei, China, 430079 Website: <http://www.servicebio.com> Tel: 027-5111-3188

Proteinase K	Servicebio	G1205
DAPI	Servicebio	G1012
anti-fade mounting medium	Servicebio	G1401
Tunel assay kit	Servicebio	G1501

2. Procedure

2.1 Deparaffinize and rehydrate: incubate sections in 2 changes of xylene, 15-20 min each. Dehydrate in 2 changes of pure ethanol for 10 min each, followed by dehydrate in gradient ethanol of 95%, 90%, 80%, and 70% ethanol, respectively, 5 min each (extend deparaffinize time slightly in winter). Wash in distilled water.

2.2 Antigen retrieval: eliminate obvious liquid, mark the objective tissue with liquid blocker pen. Add proteinase K working solution to cover objectives and incubate at 37°C for 25 min. then wash three times with PBS (pH 7.4) in a Rocker device, 5 min each. (Method for configuring working solution of proteinase K, stock solution: PBS=1:9)

2.3 Permeabilization: eliminate excess liquid, add permeabilize working solution to cover objective tissue, then incubate at room temperature for 20 min. wash three times with PBS (pH 7.4) in a Rocker device, 5 min each. (The membrane breaking fluid is 0.1% triton. Configuration method, triron stock solution: PBS=1:1000)

2.4 Equilibrium at room temperature: After the slices are slightly dried, buffer is added to the tissues in the circle, and the buffer is incubated at room temperature for 10 minutes

2.5 Tunel reaction: Take appropriate amount of TDT enzyme, dUTP and buffer in the tunel kit according to the number of slices and tissue size and mix at 1:5:50 ratio. Prepare this reaction solution according to demand before use. Add this mixture to objective tissue placed in a flat wet box, incubate at 37°C for 2 h. be sure to keep the wet box moist by adding water

2.6 DAPI counterstain in nucleus: wash three times with PBS (pH 7.4) in a Rocker device, 5 min each. Then incubate with DAPI solution at room temperature for 10 min, kept in dark place.

2.7 Mount: wash three times with PBS (pH 7.4) in a Rocker device, 5 min each. Throw away liquid slightly, then coverslip with anti-fade mounting medium.

2.8 Microscopic examination and collecting images through fluorescence microscope. DAPI emits blue light at an ultraviolet excitation wavelength of 330-380 nm and an emission wavelength of

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Wuhan Servicebio Technology Co., Ltd.

420 nm; FITC has an excitation wavelength of 465-495 nm and an emission wavelength of 515-555 nm, and emits green.

3.Results

Nucleus is blue by labeling with DAPI. TUNEL assay kit is labelled with FITC. Positive apoptosis cells are green.

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