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Add: 5th Floor, C22 Building, Biopark, No. 388 Gaoxin 2nd Road, East Lake High-tech Developing Zone, Wuhan, Hubei, China, 430079 Website: <http://www.servicebio.com> Tel: 027-5111-3188

Permeabilize solution	Servicebio	G1204
Proteinase K	Servicebio	G1205
3% H ₂ O ₂	Servicebio	G0115
Hematoxylin staining solution	Servicebio	G1004
Hematoxylin differentiate solution	Servicebio	G1039
Hematoxylin bluing solution	Servicebio	G1040
Resin mounting medium	SCRC	10004160
Tunel Assay Kit	Servicebio	G1507
DAB reagent	Servicebio	G1212

2. Procedure

2.1 Deparaffinize and rehydrate: incubate sections in 2 changes of xylene, 20 min each. Dehydrate in 2 changes of pure ethanol for 10 min each, followed by dehydrate in gradient ethanol of 95%, 90%, 80%, and 70% ethanol, respectively, 5 min each (extend deparaffinize time slightly in winter). Wash in distilled water.

2.2 Antigen retrieval: eliminate obvious liquid, mark the objective tissue with liquid blocker pen. Add proteinase K working solution to cover objectives and incubate at 37°C for 20 min. then wash three times with PBS (pH 7.4) in a Rocker device, 5 min each.

2.3 (Optional step) Permeabilization: eliminate excess liquid, add permeabilize working solution to cover objective tissue, then incubate at room temperature for 20 min. wash three times with PBS (pH 7.4) in a Rocker device, 5 min each.

2.4 Block endogenous peroxidase: wash three times with PBS (pH 7.4) in a Rocker device, 5 min each. Immerse in 3% H₂O₂ and incubate at room temperature for 20 min, kept in dark place. Then wash again three times with PBS (pH 7.4) in a Rocker device, 5 min each.

2.5 Equilibrium at room temperature: After the slices are slightly dried, buffer is added to the tissues in the circle, and the buffer is incubated at room temperature for 10 minutes

2.6 Tunel reaction: Take appropriate amount of TDT enzyme, dUTP and buffer in the tunel kit according to the number of slices and tissue size and mix at 1:5:50 ratio. Prepare this reaction solution according to demand before use. Add this mixture to objective tissue placed in a flat wet box, incubate at 37°C for 1 h. be sure to keep the wet box moist by adding water.

2.7 Reagent Streptavidin-HRP and TBST are mixed at a ratio of 1:200, added to the circle to cover the tissue, the slices are placed in a wet box, and incubated in a 37°C incubator for 30 minutes.

The slides were placed in PBS (PH7.4) and washed with shaking on a decolorizing shaker 3 times, 5 min each time.

2.8 DAB developing: dry sections slightly, add fresh prepared DAB chromogenic reagent to

marked tissue. Manage reaction time by observing in microscopy until nucleus shows brown-yellow. Then stop developing reaction by wash in pure water.

2.9 Counterstain in nucleus with Hematoxylin staining solution for 1 min and wash in pure water. Treat with the differentiate solution for a few seconds, wash in pure water. Back to blue by bluing solution, wash in pure water.

2.10 Dehydrate the slices in 4 cylinders of 100% alcohol for 5 minutes each, then soak them in n-butanol for 5 minutes, Put it in xylene for transparency, more than 5min. Dry briefly and mount with resin mounting medium.

3. Results

Nucleus stained with hematoxylin are blue. The positive apoptosis cells developed by DAB reagent have brown-yellow nucleus.