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Add: 5th Floor, C22 Building, Biopark, No. 388 Gaoxin 2nd Road, East Lake High-tech Developing Zone, Wuhan, Hubei, China, 430079 Website: <http://www.servicebio.com> Tel: 027-5111-3188

Tunel & Immunofluorescence Double Staining Protocol (Paraffin Sections)

1. Apparatus and reagents

1.1 Apparatus

Name	Producer	Model
Dehydrator	DIAPATH	Donatello
Paraffin embedding machine	武汉俊杰电子有限公司	JB-P5
pathology slicer	Leica	RM2016
Frozen machine	武汉俊杰电子有限公司	JB-L5
Water bath-slide drier	浙江省金华市科迪仪器设备有限公司	KD-P
Oven	上海慧泰仪器制造有限公司	DHG-9140A
Glass microscope slides	Servicebio	
Coverslips	CITOTEST	10212432C
Microwave	Glanze	P70D20TL-P4
Rocker	Servicebio	TSY-B
Vortex	Servicebio	MX-F
Micro-centrifuge	Servicebio	D1008E
Pipettor	Dragon	KE0003087/KA0056573
Liquid blocker pen	Servicebio	G6100
Ortho-Fluorescent Microscopy	Nikon	NIKON ECLIPSE C1
Imaging system	Nikon	NIKON DS-U3

1.2 Major reagents

Name	Producer	Code
Ethanol	SCRC	100092183
Biodewax and Clear Solution	SERVICEBIO	G1128
PBS solution	Servicebio	G0002

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Add: 5th Floor, C22 Building, Biopark, No. 388 Gaoxin 2nd Road, East Lake High-tech Developing Zone, Wuhan, Hubei, China, 430079 Website: <http://www.servicebio.com> Tel: 027-5111-3188

Permeabilize solution	Servicebio	G1204
Proteinase K	Servicebio	G1205
BSA	Servicebio	G5001
Tunel Assay Kit	Servicebio	G1501
Primary antibody		
Fluorescent-labelled secondary antibody		
DAPI	Servicebio	G1012
anti-fade mounting medium	Servicebio	G1401

2. Procedure

2.1 Deparaffinize and rehydrate: incubate sections in 2 changes of xylene, 15-20 min each. Dehydrate in 2 changes of pure ethanol for 10 min each, followed by dehydrate in gradient ethanol of 95%, 90%, 80%, and 70% ethanol, respectively, 5 min each (extend deparaffinize time slightly in winter). Wash in distilled water.

2.2 Antigen retrieval: eliminate obvious liquid, mark the objective tissue with liquid blocker pen. Add proteinase K working solution to cover objectives and incubate at 37°C for 25 min. then wash three times with PBS (pH 7.4) in a Rocker device, 5 min each. (Method for configuring working solution of proteinase K, stock solution: PBS=1:9)

2.3 Permeabilization: eliminate excess liquid, add permeabilize working solution to cover objective tissue, then incubate at room temperature for 20 min. wash three times with PBS (pH 7.4) in a Rocker device, 5 min each. (The membrane breaking fluid is 0.1% triton. Configuration method, triron stock solution: PBS=1:1000)

2.4 Equilibrium at room temperature: After the slices are slightly dried, buffer is added to the tissues in the circle, and the buffer is incubated at room temperature for 10 minutes

2.5 Tunel reaction: Take appropriate amount of TDT enzyme, dUTP and buffer in the tunel kit according to the number of slices and tissue size and mix at 1:5:50 ratio。 Prepare this reaction solution according to demand before use. Add this mixture to objective tissue placed in a flat wet box, incubate at 37°C for 2 h. be sure to keep the wet box moist by adding water.

2.6 BSA blocking: wash slides three times with PBS (pH 7.4) in a Rocker device, 5 min each. Add 3% BSA to cover the marked tissue to block non-specific binding for 30 min BSA.

2.7 Primary antibody: throw away the blocking solution slightly. Incubate slides with primary antibody (diluted with PBS appropriately) overnight at 4 °C, placed in a wet box containing a little

water.

2.8 Secondary antibody: wash slides three times with PBS (pH 7.4) in a Rocker device, 5 min each. Then throw away liquid slightly. Cover objective tissue with secondary antibody (appropriately respond to primary antibody in species), incubate at room temperature for 50 min in dark condition.

2.9 DAPI counterstain in nucleus: wash three times with PBS (pH 7.4) in a Rocker device, 5 min each. Then incubate with DAPI solution at room temperature for 10 min, kept in dark place.

2.10 Mount: wash three times with PBS (pH 7.4) in a Rocker device, 5 min each. Throw away liquid slightly, then coverslip with anti-fade mounting medium.

2.11 Microscopy detection and collect images by Fluorescent Microscopy. DAPI glows blue by UV excitation wavelength 330-380 nm and emission wavelength 420 nm; FITC glows green by excitation wavelength 465-495 nm and emission wavelength 515-555 nm; CY3 glows red by excitation wavelength 510-560 nm and emission wavelength 590 nm.

3. Results

The nuclei stained by DAPI are blue under the excitation of ultraviolet light, the tunel kit is labeled with FITC fluorescein, the positive apoptotic cells are labeled with green light, and the primary antibody is labeled with red.

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Wuhan Servicebio Technology Co., Ltd.