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Add: 5th Floor, C22 Building, Biopark, No. 388 Gaoxin 2nd Road, East Lake High-tech Developing Zone,
Wuhan, Hubei, China, 430079 Website: <http://www.servicebio.com> Tel: 027-5111-3188

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Tunel & Immunofluorescence Double Staining Protocol (Frozen-Slides)

1. Apparatus and reagents

1.1 Apparatus

Name	Producer	Model
Freezing microtome	Thermo	Cryotome E
Glass microscope slides	Servicebio	
Coverslips	CITOTEST	10212432C
Rocker	Servicebio	TSY-B
Vortex	Servicebio	MX-F
Micro-centrifuge	Servicebio	D1008E
Pipettor	Dragon	KE0003087/KA0056573
Liquid blocker pen	Servicebio	G6100
Refrigerator	Haier	BCD-192TGN
Ortho-Fluorescent Microscopy	Nikon	NIKON ECLIPSE C1
Imaging system	Nikon	NIKON DS-U3

1.2 Major reagents

Name	Producer	Model
4% paraformaldehyde	Servicebio	G1101
PBS solution	Servicebio	G0002
Proteinase K	Servicebio	G1205
Permeabilize solution	Servicebio	G1204
BSA	Servicebio	G5001
Tunel assay kit	Servicebio	G1501
Primary antibody		
Fluorescent-labelled secondary antibody		
DAPI	Servicebio	G1012
anti-fade mounting medium	Servicebio	G1401

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Add: 5th Floor, C22 Building, Biopark, No. 388 Gaoxin 2nd Road, East Lake High-tech Developing Zone, Wuhan, Hubei, China, 430079 Website: <http://www.servicebio.com> Tel: 027-5111-3188

2. Procedure

2.1 Fix frozen-slides: restore frozen-slides to room temperature and dry in air. Fix in cold acetone 10 min and then dry in air. Wash three times with PBS (pH 7.4) in a Rocker device, 5 min each.

2.2 Antigen retrieval: eliminate obvious liquid, mark the objective tissue with liquid blocker pen. Add proteinase K working solution to cover objectives and incubate at 37°C for 25 min. then wash three times with PBS (pH 7.4) in a Rocker device, 5 min each. (Method for configuring working solution of proteinase K, stock solution: PBS=1:9)

2.3 Permeabilization: eliminate excess liquid, add permeabilize working solution to cover objective tissue, then incubate at room temperature for 20 min. wash three times with PBS (pH 7.4) in a Rocker device, 5 min each. (The membrane breaking fluid is 0.1% triton. Configuration method, triton stock solution: PBS=1:1000)

2.4 Equilibrium at room temperature: After the slices are slightly dried, buffer is added to the tissues in the circle, and the buffer is incubated at room temperature for 10 minutes

2.5 Tunel reaction: Take appropriate amount of TDT enzyme, dUTP and buffer in the tunel kit according to the number of slices and tissue size and mix at 1:5:50 ratio. Prepare this reaction solution according to demand before use. Add this mixture to objective tissue placed in a flat wet box, incubate at 37°C for 2 h. be sure to keep the wet box moist by adding water.

2.6 BSA blocking: wash slides three times with PBS (pH 7.4) in a Rocker device, 5 min each. Add 3% BSA to cover the marked tissue to block non-specific binding for 30 min BSA.

2.7 Primary antibody: throw away the blocking solution slightly. Incubate slides with primary antibody (diluted with PBS appropriately) overnight at 4 °C, placed in a wet box containing a little water.

2.8 Secondary antibody: wash slides three times with PBS (pH 7.4) in a Rocker device, 5 min each. Then throw away liquid slightly. Cover objective tissue with secondary antibody (appropriately respond to primary antibody in species), incubate at room temperature for 50 min in dark condition.

2.9 DAPI counterstain in nucleus: wash three times with PBS (pH 7.4) in a Rocker device, 5 min each. Then incubate with DAPI solution at room temperature for 10 min, kept in dark place.

2.10 Mount: wash three times with PBS (pH 7.4) in a Rocker device, 5 min each. Throw away liquid slightly, then coverslip with anti-fade mounting medium.

2.11 Microscopy detection and collect images by Fluorescent Microscopy. DAPI glows blue by UV excitation wavelength 330-380 nm and emission wavelength 420 nm; FITC glows green by excitation wavelength 465-495 nm and emission wavelength 515-555 nm; CY3 glows red by excitation wavelength 510-560 nm and emission wavelength 590 nm.

3. Results

The nuclei stained by DAPI are blue under the excitation of ultraviolet light, the tunel kit is

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Wuhan Servicebio Technology Co., Ltd.

labeled with FITC fluorescein, the positive apoptotic cells are labeled with green light, and the primary antibody is labeled with red.

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