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Add: 5th Floor, C22 Building, Biopark, No. 388 Gaoxin 2nd Road, East Lake High-tech Developing Zone,
Wuhan, Hubei, China, 430079 Website: <http://www.servicebio.com> Tel: 027-5111-3188

蛋白酶 K	Servicebio	G1234
苏木素染液	Servicebio	G1004
超净快干封片胶	Servicebio	G1404
杂交缓冲液	Servicebio	G3016-3
anti-DIG-HRP	jackson	200-032-156
DAB 显色剂	Servicebio	G1211

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1. Apparatus and Reagents

1.1 Apparatus

Name	Producer	Model
Frozen slicer	Thermo	Cryotome E
Glass slide	Servicebio	
Coverslip	CITOTEST	10212432C
Enzyme-free centrifuge tubes	Servicebio	EP-150-M
Shaker	Servicebio	TSY-B
Vortex	Servicebio	MV-100
Pipettor	Dragon	KE0003087/KA0056573
Liquid Blocker PAP Pen	Servicebio	G6100
Refrigerator	HAIER	BCD-192TGN
Microscopy	Nikon	NIKON ECLIPSE CI
Imaging system	Nikon	NIKON DS-U3
Incubator	LABOTERY	GSP-70
autoclave	PANASONIC	MLS-3751L-PC

1.2 Major reagents

Reagent	Manufacturer	Article Number
DEPC	Amresco	E174
4% of paraformaldehyde (DEPC water)	Servicebio	G1113
OCT embedding agent	Servicebio	G6059-110ML
Ethanol	SCRC	100092683
BioDewax and Clear solution	Servicebio	G1128
Methanol	SCRC	40030561
Hydrogen peroxide	SCRC	10011208
Sucrose	Servicebio	G5031
PBS	Servicebio	G0002
20×SSC solution	Servicebio	G3016-4
BSA	Servicebio	G5001

Proteinase K	Servicebio	G1234
Hematoxylin dye	Servicebio	G1004
SweSuper Clean BioMount Medium	Servicebio	G1404
Hybridization buffer	Servicebio	G3016-3
anti-DIG-HRP	Jackson	G3016-3
DAB reagent	Servicebio	G1211

2. The Steps of the Experiment

2.1. **Tissue fixation:** take out the organization, wash clean, then Immediately put in the fixed fluid (DEPC) above 12h.

2.2. **Tissue dehydration:** Transfer tissue to 15% sucrose solution for 8 h and were subsequently to a 30% sucrose solution overnight.

2.3. **Frozen section fixation:** The frozen sections were dried at room temperature, fixed in 4% paraformaldehyde (DEPC) for 10 min, and shaken on the shaker in PBS (pH 7.4) for three times, 5 min each.

2.4. **Digestion:** Mark the objective tissue with liquid blocker pen, according to the characteristics of different tissues and different indicators, add proteinase K (20 ug/ml) to cover tissues and incubate at 37°C for _____minutes. Washed with sterilized water and then washed three times in PBS for 5 min each time.

2.5. **Block endogenous peroxidase:** add 3% methanol-H₂O₂, incubate in dark at room temperature for 15min. Wash slides in PBS(PH7.4) three times for 5 min each with gentle agitation.

2.6. **Prehybridization:** add hybridization buffer onto specimen and incubate at 37°C for 1h.

2.7. **Hybridization:** remove the pre-hybridization solution, add the _____ probe hybridization solution with concentration of _____, and incubate the section in a humidity chamber and hybridize overnight at _____°C.

2.8. **Washing:** Remove the hybridization solution. Wash sections in 2×SSC for 10 min at 37°C , 1×SSC two times for 5 min each at 37°C, and 0.5×SSC for 10 min at room temperature. Formamide washing can be added if there are more non-specific hybrids.

2.9. **Blocking:** add blocking solution (Rabbit serum) to the section and incubate at room temperature for 30 min.

2.10. Add the mouse anti-digoxigenin-labeled peroxidase (anti-DIG-HRP): Remove the blocking solution and add anti-DIG-HRP. Incubate at 37 °C for 50 min, then Wash sections in PBS three times for 5 min each.

2.11. **DAB developing:** dry sections slightly, add fresh prepared DAB chromogenic reagent to marked tissue. Manage reaction time by observing in microscopy until positive expression appears

brown-yellow. Then stop developing reaction by wash in running tap water.

2.12. **Stain cell nuclei (counter stain):** Counterstain in nucleus with hematoxylin staining solution for 3 min and wash in tap water. Treat with the differentiate solution for a few seconds, wash in running tap water. Back to blue by bluing solution, wash in running tap water.

2.13. **Mounting:** Dehydrate successively in gradient ethanol of 75%, 85%, and 2 changes of pure ethanol, respectively, 6 min each. Clear in BioDewax and Clear solution for 6 min and mount with SweSuper Clean BioMount Medium.

2.14. Image the results using a brightfield fluorescence microscope. Store slides at RT.

3. Interpretation of the results

The nucleus stained with hematoxylin are blue, and the positive expression of DAB was brownish yellow. The results of in situ hybridization of mRNA were cytoplasmic positive and a few nuclear positive were normal. MicRNA and lncRNA were expressed differently. According to the expression of color depth.

Note: All reagents, instrument need RNase free.

Attached table 1 probe information.

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Wuhan Servicebio Technology Co., Ltd.