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(Red) Tunel Assay Protocol (Frozen-Slides)

1. Apparatus and reagents

1.1 Apparatus

Name	Producer	Model
Freezing microtome	Thermo	Cryotome E
Glass microscope slides	Servicebio	
Coverslips	CITOTEST	10212432C
Rocker	Servicebio	TSY-B
Vortex	Servicebio	MX-F
Micro-centrifuge	Servicebio	D1008E
Pipettor	Dragon	KE0003087/KA0056573
Liquid blocker pen	Servicebio	G6100
Refrigerator	Haier	BCD-192TGN
Ortho-Fluorescent Microscopy	Nikon	NIKON ECLIPSE C1
Imaging system	Nikon	NIKON DS-U3

1.2 Major reagents

Name	Producer	Model
4% paraformaldehyde	Servicebio	G1101
PBS solution	Servicebio	G0002
Proteinase K	Servicebio	G1205
Permeabilize solution	Servicebio	G1204
DAPI	Servicebio	G1012
anti-fade mounting medium	Servicebio	G1401
Tunel assay kit	Servicebio	G1502

2. Procedure

2.1 Fix frozen-slides: restore frozen-slides to room temperature and dry in air. Fix in cold acetone 10 min and then dry in air. Wash three times with PBS (pH 7.4) in a Rocker device, 5 min each

2.2 Antigen retrieval: eliminate obvious liquid, mark the objective tissue with liquid blocker pen. Add proteinase K working solution to cover objectives and incubate at 37°C for 25 min. then wash three times with PBS (pH 7.4) in a Rocker device, 5 min each. (Method for configuring working solution of proteinase K, stock solution: PBS=1:9)

2.3 Permeabilization: eliminate excess liquid, add permeabilize working solution to cover objective tissue, then incubate at room temperature for 20 min. wash three times with PBS (pH 7.4) in a Rocker device, 5 min each. (The membrane breaking fluid is 0.1% triton. Configuration method, triton stock solution: PBS=1:1000)

2.4 Equilibrium at room temperature: After the slices are slightly dried, buffer is added to the tissues in the circle, and the buffer is incubated at room temperature for 10 minutes

2.5 Tunel reaction: Take appropriate amount of TDT enzyme, dUTP and buffer in the tunel kit according to the number of slices and tissue size and mix at 1:5:50 ratio. Prepare this reaction solution according to demand before use. Add this mixture to objective tissue placed in a flat wet box, incubate at 37°C for 2 h. be sure to keep the wet box moist by adding water

2.6 DAPI counterstain in nucleus: wash three times with PBS (pH 7.4) in a Rocker device, 5 min each. Then incubate with DAPI solution at room temperature for 10 min, kept in dark place.

2.7 Mount: wash three times with PBS (pH 7.4) in a Rocker device, 5 min each. Throw away liquid slightly, then coverslip with anti-fade mounting medium.

2.8 Microscopy detection and collect images by Fluorescent Microscopy. DAPI glows blue by UV excitation wavelength 330-380 nm and emission wavelength 420 nm; CY3 glows red by excitation wavelength 510-560 nm and emission wavelength 590 nm by excitation wavelength 465-495 nm and emission wavelength 515-555 nm;

3. Results

The nuclei stained by DAPI are blue under the excitation of ultraviolet light, the tune3 kit is labeled with CY3 fluorescein, and the positive apoptotic nuclei are red