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<http://www.servicebio.cn>

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Add: 5th Floor, C22 Building, Biopark, No. 388 Gaoxin 2nd Road, East Lake High-tech Developing Zone, Wuhan, Hubei, China, 430079 Website: <http://www.servicebio.com> Tel: 027-5111-3188

Experimental Report of Caspary Band Staining

1. Apparatus and Reagents

1.1 Major apparatus

Brand name	Manufactures	Model
Glass slide	Servicebio	GFL-230
Cover glass	Jiangsu Shitai Experimental Equipment Co.,Ltd	10127105P-G
Upright optical microscope	Nikon	Nikon Eclipse E100
Imaging system	Nikon	Nikon DS-U3
Dehydrator	DIAPATH	Donatello
Embedding machine	Wuhan Junjie Electronics Co., Ltd	JB-P5
Pathology slicer	Leica	RM2016
Frozen platform	Wuhan Junjie Electronics Co., Ltd	JB-L5
Organizer	KEDEE	KD-P
oven	Labotery	GFL-230

1.2 Major reagents

Reagent	Manufactures	Cat
Aniline blue staining solution	Servicebio	71003644
Berberine hydrochloride	Macklin	B802465
Neutral resin	Sinopharm Chemical Reagent Co.,Ltd	10004160

2. Experimental

2.1 Dewaxing as followed:

Xylene I for 20 min;
Xylene II for 20 min;
100% ethanol I for 5 min;
100% ethanol II for 5 min;
75% ethanol for 5 min;
Rinsing with tap water.

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Wuhan Servicebio Technology Co., Ltd.

2.2 Put the slides into 0.1% Berberine hydrochloride solution for 1 hour under dark conditions. Wash the slides with detergent for 3 times. Put the slides into 0.05% aniline blue solution for 30 min under dark conditions. Wash the slides with detergent for 3 to 5 times and the filter paper absorbs the excess liquid.

2.3 Seal: Drop anti-fluorescence quenching sealing film liquid sealing sheet, cover glass slide.

2.4 Take photos: observe and take photos. Observe the deposition of Caspary band with fluorescence microscope (DAPI filter). The fluorescence is easy to quench. Take photos as soon as possible and avoid light. The slices can be stored at 4 °C for a short time.

3. The results were as follows

On the upper, lower, and radial walls of the inner cortex cells, bright blue ribbons of The Caspary band encircle the cells.

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