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Scanning Electron Microscopy-Energy Dispersive X-Ray Spectroscopy Report

1. Apparatus and Reagents

1.1 Major Apparatus

Name	Producer	Model
Critical Point Dryer	Quorum	K850
Lon Sputtering Apparatus	HITACHI	MC1000
Scanning Electron Microscope	HITACHI	SU8100
EDS	Oxford	ULTIMMAX 100

1.2 Major Reagents

Name	Producer	Code
Fixative for TEM	Servicebio	G1102
Ethanol	Sinaopharm Group Chemical Reagent Co. LTD	100092183
Isoamyl acetate	Sinaopharm Group Chemical Reagent Co. LTD	10003128
PBS	Servicebio	G0002
OsO ₄	Ted Pella Inc	

2. Procedure

2.1 Harvest tissue block and fixation: Targeted fresh tissues should be selected to minimize mechanical damage such as pulling, contusion and extrusion. Use a sharp blade to cut and harvest fresh tissue blocks quickly within 1-3 minutes. The area of tissue block should be no more than 3 mm². Wash tissues with PBS gently to remove the blood and hair, etc. Label the target side of tissue (the side you want to observe) by any way, such as making cuts on the opposite side. Make sure to protect tissue blocks, especially the target side, from mechanical damage such as forceps extrusion. The washed tissue blocks are immediately fixed by electron microscopy fixative for 2 hours at room temperature, then transferred into 4°C for preservation and transportation.

For adherent cell: Seed cells on a sterile cover glass in a petri dish. Remove the culture medium, then wash slide gently with PBS, followed by adding electron microscopy fixative into petri dish. After fixing for 2 hours at room temperature, transfer the petri dish to 4°C for preservation and transport. Note to avoid sever shock which may result in the cells dropping off from cover glass.

2.2 Post-fix: Wash tissue blocks with 0.1 M PB (pH 7.4) for 3 times, 15 min each. Then transfer tissue blocks into 1% OsO₄ in 0.1 M PB (pH 7.4) for 1-2 h at room temperature. After that, wash tissue blocks in 0.1M PB (pH 7.4) for 3 times, 15 min each.

2.3 Dehydrate as followed:

30% ethanol for 15 min;
50% ethanol for 15 min;
70% ethanol for 15 min;
80% ethanol for 15 min;
90% ethanol for 15 min;
95% ethanol for 15 min;
Two changes of 100% ethanol for 15 min;
Finally, isoamyl acetate for 15 min.

2.4 Drying: Dry samples with Critical Point Dryer.

Note: Dry samples and inorganic materials should ignore all the steps above and do the following step directly.

2.5 Conductive metal coating: Specimens are attached to metallic stubs using carbon stickers and sputter-coated with gold for 30s.

2.6 Elemental analysis : The element analysis of the area to be measured was observed under scanning electron microscope.

Point scanning: Three points were randomly selected from the area to be tested, energy signals were collected, element content of each point was detected, and atlas and percentage table were saved.

Line scan: A linear trajectory was randomly selected in the region to be tested, and the change of element concentration on the line was detected. The fluctuation of the curve reflects the change in the concentration of the corresponding element along the scan line.

Surface scan: A certain area on the surface of the sample is used as the detection area. By collecting characteristic X-ray photons, the plane distribution diagram of each element can be obtained, and each element is represented by different colors.

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Wuhan Servicebio Technology Co., Ltd.