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1.2. Major reagents

reagent	manufacturer	article number
Ethanol	SCRC	100092683
BioDewax and Clear solution	Servicebio	G1128
PBS	Servicebio	G0002
Proteinase K	Servicebio	G1234
4% of paraformaldehyde (DEPC water)	Servicebio	G1113
Methanol	SCRC	
Hydrogen peroxide	SCRC	
DEPC	Amresco	E174
20×SSC solution	Servicebio	G3016-4
BSA	Servicebio	G5001
Hematoxylin dye	Servicebio	G1004
SweSuper Clean BioMount Medium	Servicebio	G1404

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Add: 5th Floor, C22 Building, Biopark, No. 388 Gaoxin 2nd Road, East Lake High-tech Developing Zone, Wuhan, Hubei, China, 430079 Website: <http://www.servicebio.com> Tel: 027-5111-3188

hybridization buffer	Servicebio	G3016-3
anti-DIG-HRP	Jackson	
DAB reagent	Servicebio	G1211

2. The steps of the experiment

2.1. **Organization fixed:** Take out the organization, Wash clean, then Immediately put in the fixed fluid (DEPC) to fix 2-12h.

2.2. **Dehydration:** The tissue is dehydrated by gradient alcohol, paraffin, Embedding.

2.3. **Section:** The paraffin is sliced through the slicer, the piece of the slice machine and the 62 - degree oven roast for 2 hours.

2.4. **Dewaxing and dehydration:** Soak sections in 2 changes of Dewaxing Transparent Liquid, 15 minutes each. Dehydrate in 2 changes of pure ethanol for 5 minutes each. Then, followed respectively by dehydrating in gradient ethanol of 85% and 75% ethanol 5 minutes each. Wash in DEPC dilution.

2.5. **Digestion:** According to the tissue fixation time, the slices are boiled in the retrieval solution for 10-15 minutes and naturally cooled. Mark the objective tissue with liquid blocker pen, according to the characteristics of tissues, Add proteinase K (20ug/ml) working solution to cover objectives and incubate at 37°C for ___ min. Washing with pure water, then wash three times with PBS (pH 7.4) in a Rocker device, 5 min each.

2.6. **Block endogenous peroxidase:** add 3% methanol-H₂O₂ , incubate in dark at room temperature for 15min. Wash slides in PBS (pH 7.4) three times for 3 min each, with gentle agitation and protected from light.

2.7. **Pre-hybridization:** Add Pre-hybridization solution to each section and incubate for 1 h at 37°C.

2.8. **Hybridization:** Discard the pre-hybridization solution, add the probe hybridization solution , concentration ___, and incubate the section in a humidity chamber and hybridize overnight at ___ °C.

2.9. **Washing:** Remove the hybridization solution. Wash sections with 2×SSC for 5 min at 37°C , Wash sections in 1×SSC two times for 5 min each at 37°C, and wash in 0.5×SSC for 10 min at room temperature. Formamide washing can be added if there are more non-specific hybrids.

2.10. **Imaging oligo (DIG):** Add the hybridization solution containing Imaging Oligo (DIG), Dilution ratio 1:400. 42 ° FOR 3 hours. Wash sections with 2×SSC for 5 min at 37°C, Wash sections in 1×SSC two times for 5 min each at 37°C, and wash in 0.5×SSC for 10 min at room

temperature.

2.11. **Blocking:** add blocking solution (Rabbit serum) to the section and incubate at room temperature for 30 min.

2.12. Add the mouse anti-digoxigenin-labeled peroxidase (anti-DIG-HRP): Remove the blocking solution and add anti-DIG-HRP. Incubate at 37 °C for 40 min, then Wash sections in PBS four times for 5 min each.

2.13. **DAB developing:** dry sections slightly, add fresh prepared DAB chromogenic reagent to marked tissue. Manage reaction time by observing in microscopy until positive expression appears brown-yellow. Then stop developing reaction by wash in running tap water.

2.14. Counterstain in nucleus with Hematoxylin staining solution for 3 min and wash in tap water. Treat with the differentiate solution for a few seconds, wash in running tap water. Back to blue by bluing solution, wash in running tap water.

2.15. Dehydrate successively in gradient ethanol of 75%, 85%, and 2 changes of pure ethanol, respectively, 6 min each. Clear in Dewaxing Transparent Liquid for 6 min and mount with SweSuper Clean BioMount Medium.

2.16. Image the results using a brightfield fluorescence microscope. Store slides at RT.

3. Interpretation of the results

The nucleus stained with hematoxylin are blue, and the positive expression of DAB was brownish yellow. The results of in situ hybridization of mRNA were cytoplasmic positive and a few nuclear positive were normal. MicRNA and lncRNA were expressed differently. According to the expression of color depth.

Note: All reagents, instrument need RNase free.

Attached table 1 probe information.