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Victoria Blue Staining Experimental Report

1. Lab Equipment and Reagents

1.1 Lab equipment

Items	Manufacturer	Model
Dehydrator	DIAPATH	Donatello
embedding machine	Wuhan Junjie Electronics Co., Ltd.	JB-P5
Pathology microtome	Shanghai Leica Instruments Co., Ltd.	RM2016
Frozen platform	Wuhan Junjie Electronics Co., Ltd.	JB-L5
Water Bath-Slide Drier	Zhejiang Jinhua Kedi Instrumental Equipment CO.,LTD	KD-P
Laboratory oven	Tianjin Labotery Instrument Equipment Co., Ltd.	GFL-230
Microscope slide &cover glass	Citotest Labware Manufacturing Co.,Ltd	10127105P-G
Upright optical microscope	Nikon Japan	Nikon Eclipse E100
Imaging system	Nikon Japan	NIKON DS-U3

1.2 Chemical Reagents

Items	Manufacturer	Model
Absolute alcohol	Sinopharm Chemical Reagent Co., Ltd.	100092683
Xylene	Sinopharm Chemical Reagent Co., Ltd.	10023418
Victoria blue staining kit	Servicebio	G1055
Neutral balsam	Sinopharm Chemical Reagent Co., Ltd.	10004160

2. Experimental steps

2.1 Paraffin section deparaffinization and rehydration: put the slides into xylene I 20minutes-xylene II 20 minutes-absolute ethanol I 5 min-absolute ethanol II 5 min-75% alcohol for 5 min, tap water washing.

2.2 Tissue acidification: draw a circle to enclose the tissue, add Victoria blue acidification solution (Victoria Blue A and Victoria Blue B 1:1 and mix, ready for usage), oxidize the tissue for 5 minutes, dip the slides in two cylinders of water for 10s, and dry the slides to remove excess water.

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Add: 5th Floor, C22 Building, Biopark, No. 388 Gaoxin 2nd Road, East Lake High-tech Developing Zone, Wuhan, Hubei, China, 430079 Website: <http://www.servicebio.com> Tel: 027-5111-3188

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Wuhan Servicebio Technology Co., Ltd.

2.3 Tissue bleaching: The tissues are bleached with Victoria Blue C for 2 minutes, put the slides soak in 3 cylinders of water for 5 s, and dry the slides to remove excess water.

2.4 Tissue coloring: dip into 70% ethanol, put the slides into Victoria blue staining D for 24h (covered), dip into 70% ethanol washing for 2 times, 10s each time, until there is no staining solution attached to the slide surface, running water washing.

2.5 Tissue differentiation: checking the degree of staining of elastic fibers under a microscope, if deep coloring, put slides into 75% ethanol for 5s, tap water washing to terminate differentiation microscope examination, repeatedly differentiation, washing and microscope examination, until elastic fiber is blue, light blue background or nearly colorless;

2.6 Counterstaining: staining with Victoria Blue Staining Solution E for 1-5 minutes and running water washing.

2.7 Dehydration and sealing: put the slide into absolute ethanol I 5min-absolute ethanol II 5min-absolute ethanol III 5min-xylene I 5min-xylene II 5min transparent, neutral balsam sealing.

2.8 Microscope examination, image collection and analysis.

3. Results

The elastic fiber is blue, nucleus is red.