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1. Apparatus and Reagents

1.1 Major apparatus

Name	Producer	Model
Frozen slicer	Thermo	CRYOSTAR NX50
Slicer	LEICA	LEICA 819
Glass slide	Servicebio	

1.2 Major reagents

Name	Producer	Code
Fixative	Servicebio	G1101
Oil Red solution	Servicebio	G1015
Hematoxylin solution	servicebio	G1004
Differentiation solution (With 60% ethanol as solvent)	servicebio	G1039
Scott Tap Bluing	servicebio	G1040
Glycerin gelatin sealed tablets	Servicebio	G1402
Isopropanol	SCRC	80109218

2. Procedure

2.1 Reheat and dry the frozen slices, then fix it in the fixative solution for 15 minutes, wash with tap water, and dry.

2.2 Stain sections with Oil Red solution for 8-10 min in the dark, and cover it with lid during dyeing.

2.3 Take out the slices, stay for 3s and then immerse them in two cup of 60% isopropanol for differentiation in turn, 3s and 5s respectively. The slices were immersed to 2 cup of pure water in turn for 10s each.

2.4 Take out the slices, immerse in hematoxylin for 3-5 min after 3s, and then rinse in 3 cup of pure water for 5s, 10s, and 30s in turn. Treat it with differentiation solution (60% alcohol as solvent) for 2-8s, 2 cup of distilled water for 10s each, and Scott Tap Bluing for 1s. Then lightly dip the slices in 2 cylinders of tap water for 5s and 10s in turn, and check the staining effect by microscope.

2.5 Seal the slices with glycerin gelatin.

2.6 Observe with microscope inspection, image acquisition and analysis.

3. Results

Color	Result
Orange/Red	Lipid droplets
Blue	Nuclei

4. Precautions

4.1 If it is a fresh tissue frozen section, it needs to be fixed before staining; if the tissue is fixed after freezing section, you should stain after the slices are dried.

4.2 It should be gentle during the entire operation process to avoid fat loss or displacement.

4.3 The staining results cannot be stored for a long time, and observation and photographing should be done as soon as possible after sealing.

4.4 Glycerin gelatin is solidified at room temperature, so generally stored it in an oven at 60°C. It is also not suitable to press the slide or pull off the coverslip, if there are bubbles after seal the slice. You can put it into warm water to let the coverslip fall off by itself, and then resealed to prevent fat displacement.

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Wuhan Servicebio Technology Co., Ltd.