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Plant Trypan Blue Dyeing Experiment Report

1. Experimental Equipments and Reagents

1.1 Experimental equipments

Name	Manufactor	Model
Clean bench	Su Jing An Tai	SW-CJ-1FD
Specimen bottle	Servicebio	G6008-40ML
Camera	Canon	70D

1.2 Main experimental reagents

Reagent name	factory	Article number
Absolute ethanol	Sinopharm Group Chemical Reagent Co., Ltd.	100092683
Plant trypan blue dye solution set	Servicebio	G1020

2. Experimental steps

2.1 Take the appropriate amount of plant trypan blue dye solution A and absolute ethanol at a 1:1 ratio, use tweezers to remove the plant seedlings or plant leaves from the fixing solution, immerse them in pure water and rinse 5 times back and forth, placed on the filter paper Absorb excess water and immerse in this dilution, stain at room temperature for 2~6 h in the dark.

2.2 Take out the plant seedlings or leaves with tweezers, immerse them in pure water and rinse 5 times back and forth, put them on the filter paper to absorb excess water, and immerse them in the plant trypan blue dye solution B for 3~16h.

2.3 Take out the plant seedlings or leaves with tweezers, immerse them in pure water and rinse 5 times, put them on the filter paper to absorb excess water, store them in a specimen bottle containing 20ml plant trypan blue dye solution C at room temperature, and take pictures after 30 minutes Save the result.

3. The results were as follows

The dead cells in the seedlings or leaves of the plant or the damaged parts of the cell membrane will be dark blue, and the other parts are almost colorless.