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Wuhan Servicebio Technology Co., Ltd.

PCR amplification and Agarose Gel Electrophoresis

1. Laboratory Equipment and Reagents

1.1 Laboratory equipment

Equipment	Manufacturers	Model
Centrifuge	DragonLab	D3024R
PCR instrument	Beijing Dongsheng Innovation Biotechnology Co., Ltd.	ETC811
Constant temperature oscillation shaker	Changzhou Aohua Instrument Co., Ltd	SHZ-82A
Clean bench	Suzhou Antai Air Tech Co., Ltd	SW-CJ-1FD
Electrophoresis	Servicebio	FW-600
Gel imaging system	Tanon Science & Technology Co., Ltd.	Tanon-1600R
Standard reagent type pure water meter	Qingdao Fulum Technology Co., Ltd	FBZ2001-up-p
Ultramicro spectrophotometer	Thermo	NanoDrop2000

1.2 Reagents

Reagents	Manufacturers	Order
Servicebio Genomic DNA Kit	Servicebio	G3633
2×Fast Pfu PCR Master Mix	Servicebio	G3305
2×Fast Pfu PCR Master Mix	Servicebio	G3305
Agarose	Servicebio	G5056
GN100bp DNA Ladder II	Servicebio	G3366
Centrifuge tube(1.5 mL)	Servicebio	EP-150-J
TIP	Servicebio	
Recombinant Proteinase K	Servicebio	G1234
Red Blood Cell Lysis Buffer	Servicebio	G2015
50×TAE	Servicebio	G3001
6× DNA Loading Buffer	Servicebio	G3011
SerRed	Servicebio	G3606
Water Nuclease-Free	Servicebio	G4700
Ethanol absolute	Sinopharm Group Chemical Reagent Co., Ltd.	10009218

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Add: 5th Floor, C22 Building, Biopark, No. 388 Gaoxin 2nd Road, East Lake High-tech Developing Zone,
Wuhan, Hubei, China, 430079 Website: <http://www.servicebio.com> Tel: 027-5111-3188

2. Experimental Steps

2.1 Genomic DNA Extraction Experimental steps

2.1 Samples preparation

2.1.1 For blood, please use 50-100 μ L fresh, frozen or anticoagulant adding blood. If less than 200 μ L, please make up with **buffer GA** to 200 μ L.

2.1.2 If the sample is blood from poultry, birds, amphibians, of which red blood cells have nucleolus, the amount should be reduced to 5-20 μ L and adjust the volume to 200 μ L with **buffer GA**.

2.1.3 The adherent cells should be treated to cell suspension first, then centrifuge the cells for 5 min at 1,000 rpm, then discard the flow-through and re-suspend cell pellet in 200 μ L **buffer GA**.

2.1.4 Animal tissue (2-30mg) should be treated to cell suspension first, then centrifuge the cells for 1 min at 12,000 rpm, then discard the flow-through and re-suspend cell pellet in 200 μ L **buffer GA**.

2.2 Add 20 μ L **Proteinase K Stock Solution**, mix thoroughly by vortex. If the sample is tissue: incubate at 56°C for 30-60 min until the tissue is completely lysed (mix thoroughly by vortex per 10 min).

2.3 Add 200 μ L **Buffer GB** to the sample, mix thoroughly for 6-8 times, and incubate at 70°C for 10 min to yield a homogeneous solution. Briefly centrifuge the 1.5 mL microcentrifuge tube to remove drops from the inside of the lid.

2.4 Add 200 μ L ethanol to the sample, and mix thoroughly by vortex for 15 s. A white precipitate may form on addition of ethanol. Briefly centrifuge the 1.5 mL microcentrifuge tube to remove drops from the inside of the lid.

2.5 Pipet the mixture from step 4 into the Bind DNA Mini Columns (in a 2 mL collection tube) and centrifuge at 10000g for 30 s. Discard flow-through and place the spin column into the collection tube.

2.6 Add 500 μ L **Buffer PD** to Bind DNA Mini Columns, and centrifuge at 10000g for 1 min, then discard the flow-through and place the spin column into the collection tube.

2.7 Add 600 μ L **Buffer PW** (Ensure ethanol has been added) to Bind DNA Mini Columns, and centrifuge at 10000g for 1 min. Discard the flow-through and place the spin column into the collection tube.

2.8 Repeat Step 7.

2.9 Centrifuge at 10000g for 2 min to dry the membrane completely. Note: The residual ethanol of buffer PW may have some effect on downstream application.

2.10. Place the Bind DNA Mini Column in a new clean 1.5 mL microcentrifuge tube, and pipet 50-100 μL **Elution Buffer or ddH₂O** directly to the center of the membrane. Incubate at room temperature (15-25°C) for 2-5 min, and then centrifuge for 2 min at 10000g.

2.2 PCR amplification

2.2.1 0.2mLPCR tube, 20ul pcr amplification system.

2×Fast Pfu PCR Master Mix	25 μL
Forward Primer (10 μM)	1 μL
Reverse Primer (10 μM)	1 μL
Template	2 μL
Water Nuclease-Free	Add to 50 μL

2.2.2 PCR amplification

Pre-denaturation	98°C, 2min	} 30×cycle
Denaturation	98°C, 20s	
Annealing	55°C, 20s	
Extension	72°C, 10s	
End extension	72°C, 5min	

3. Agarose gel electrophoresis:

3.1 Melt 3g agarose in 100 mL of 1X TAE in the microwave until all agarose is dissolved and there are no stringy pieces. Cool the liquid under cold running water for 10-15 seconds or allow to cool by letting it sit at RT until it is not too hot to hold. Add a very small amount of SerRed (1ul per 20mL is sufficient). You can do this simply by dipping the pipet tip into the stock solution and then swirling this into the liquid agar. SerRed is used to visualize the DNA when viewed under UV light.

3.2 Pour the gel into the gel mold held in place by the clamp, with the desired comb in place. Allow the gel to dry for about 15 to 20 minutes.

3.3 Remove comb and transfer gel in gel mold to gel box with TAE buffer, making sure that the gel is completely submerged (do not fill past max fill line).

3.4 Load the gel. Each sample needs loading buffer, which should be used at a ratio of 1: 6. Loading buffer is used as a marker for migration rates and also to give the sample density so that it sinks to the bottom of the well.

3.5 After all samples and ladder(s) have been loaded, place the lid on the gel box, make sure it is plugged into the power, and turn the power on. Run it at about 92 volts for 30-45 minutes, or until loading buffer bands are where you want them.

3.6 Check the gel when it is done using the Gel imaging system. If it's not spread out enough, run

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Wuhan Servicebio Technology Co., Ltd.

longer. Remember to turn off the power when you disconnect the lid from the box! Also, when not in use, please make sure the lid is on to minimize evaporation of buffer and to keep dust out. When you have taken a picture or saved the image, throw the gel away in the trash.