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Add: 5th Floor, C22 Building, Biopark, No. 388 Gaoxin 2nd Road, East Lake High-tech Developing Zone, Wuhan, Hubei, China, 430079 Website: <http://www.servicebio.com> Tel: 027-5111-3188

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EndoFree Maxi Plasmid Lab Report

1. Laboratory Equipment and Reagents

1.1 Laboratory equipment

Equipment	Manufacturers	Model
Centrifuge	Heal Force	Neofuge 15R
Spectrophotometer	Thermo	NANODROP 2000
Standard reagent type pure water meter	Qingdao Fulum Technology Co., Ltd	FBZ2001-up-p
Clean bench	Suzhou Antai Air Tech Co., Ltd	SW-CJ-1FD
Constant temperature oscillation shaker	Changzhou Aohua Instrument Co., Ltd	SHZ-82A

1.2 Reagents

Reagents	Manufacturers	Order
EndoFree Maxi Plasmid Kit	Tiagen Biotech(Beijing) CO.,LTD	DP117
Anhydrous ethanol	Sinopharm Group Chemical Reagent Co., Ltd.	10009218
Luria-Bertani Liquid Medium(powder)	Servicebio	G3102
Luria-Bertani Solid Medium(powder)	Servicebio	G3101
Ampicillin	Servicebio	G5053
Kanamycin	Servicebio	G5054-1G
Isopropyl alcohol	Sinopharm Group Chemical Reagent Co., Ltd.	80109218
Centrifuge tube	Servicebio	
TIP	Servicebio	

2 EndoFree Maxi Plasmid Experimental Steps

2.1 To cultivate 200 ml bacterial culture, with 37°C、220rpm、16h.

2.2 Column equilibration: place a Spin Column CP6 into a 50 ml collection tube (supplied in the kit) and add 2.5 ml Buffer BL to Spin Column CP6. Centrifuge for 2 min at 8,000 rpm ($\sim 8,228 \times g$). Discard the flow-through, and place Spin Column CP6 into the same collection tube.

2.3 Harvest 100 ml (for low copy plasmid, please harvest 200 ml) overnight cultured bacterial cells by centrifuging at 8,000 rpm ($\sim 8,228 \times g$) for 3 min at room temperature (15-25°C), and then remove all traces of supernatant.

2.4 Try to remove all traces of supernatant, use clean filter paper to absorb the fluids inside the tube wall.

2.5 Resuspend pelleted bacterial cells in 8 ml Buffer P1 (Ensure that RNase A has been added). The bacteria should be resuspended completely by vortex or pipetting up and down until no cell clumps remain.

2.6 Add 8 ml Buffer P2 and mix thoroughly by inverting the tube 6- 8 times, then incubate at room temperature for 5 min.

2.7 Add 8 ml Buffer P4, and mix immediately and thoroughly by gently inverting 6-8 times, until the whole solution become cloudy. Incubate at room temperature for 10 min. Centrifuge for 5-10 min at 8,000 rpm ($\sim 8,228 \times g$), the white material should be in the bottom of the centrifuge tube (prolong centrifugation time properly). Transfer the supernatant into a Filtration CS1 (avoid transferring large clump into the Filtration CS1, which will clog the filtration membrane). Gently insert the plunger into the Filtration CS1 and filter the cell lysate into a new 50 ml tube (not supplied in the kit).

2.8 contamination), mix completely by reverting upside and down and then transfer all solution to the Spin Column CP6 (put Spin Column CP6 into 50 ml collection tube).

2.9 Centrifuge for 2 min at 8,000 rpm ($\sim 8,228 \times g$). Discard the flow-through and place the Spin Column CP6 back into the same collection tube.

2.10 Add 10 ml Buffer PW to the Spin Column CP6 and centrifuge at 8,000 rpm ($\sim 8,228 \times g$) for 2 min. Discard the flow-through and place the Spin Column CP6 back into the same collection tube.

2.11 Repeat step 9.

2.12 Add 3ml 96-100% ethanol to the Spin Column CP6 (put the CP6 in a collection tube). Centrifuge for 2 min at 8,000 rpm ($\sim 8,228 \times g$), discard the flow-through.

2.13 Put Spin Column CP6 back to collection tube, centrifuge at 8,000 rpm ($\sim 8,228 \times g$) for 5 min for removing residual ethanol.

2.14 To elute DNA, place the Spin Column CP6 in a clean 50 ml collection tube (supplied in the kit) and add 1-2 ml Buffer TB to the center of the membrane and incubate 5 min at room temperature, centrifuge at 8,000 rpm ($\sim 8,228 \times g$) for 2 min. Transfer the eluate from 50 ml centrifuge tube to a clean 1.5 C for storage. °ml centrifuge tube and put at -20 for storage.

2.2 Determination of Plasmid Concentration: Ssee Excel.