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Add: 5th Floor, C22 Building, Biopark, No. 388 Gaoxin 2nd Road, East Lake High-tech Developing Zone, Wuhan, Hubei, China, 430079 Website: <http://www.servicebio.com> Tel: 027-5111-3188

Immunofluorescence Protocol (Paraffin-slides)

1. Apparatus and Reagents

1.1 Apparatus

Name	Producer	Model
Dehydrator	DIAPATH	Donatello
Paraffin embedding machine	Wuhan Junjie Electronics Co., Ltd.	JB-P5
pathology slicer	Leica	RM2016
Frozen machine	Wuhan Junjie Electronics Co., Ltd.	JB-L5
Water bath-slide drier	Zhejiang Jinhua Kedi Instrumental Equipment Co., Ltd.	KD-P
Oven	Shanghai Huitai Instrument Manufacturing Co., Ltd.	DHG-9140A
Glass microscope slides	Servicebio	
Coverslips	CITOTEST	10212432C
Microwave	Glanze	P70D20TL-P4
Rocker	Servicebio	TSY-B
Vortex	Servicebio	MX-F
Micro-centrifuge	Servicebio	D1008E
Pipettor	Dragon	KE0003087/KA0056573
Liquid blocker pen	Servicebio	G6100
Slice scanner	Pannoramic	Pannoramic MIDI:3Dhistech

1.2 Major reagents

Name	Producer	Code	Dilution
Ethanol	SCRC	100092183	
Biodewax and Clear Solution	SERVICEBIO	G1128	
EDTA antigen retrieval solution (pH 8.0)	Servicebio	1330-20-7 G1206	
EDTA antigen retrieval solution (pH 9.0)	Servicebio	G1203	

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Add: 5th Floor, C22 Building, Biopark, No. 388 Gaoxin 2nd Road, East Lake High-tech Developing Zone,
Wuhan, Hubei, China, 430079 Website: <http://www.servicebio.com> Tel: 027-5111-3188

Sodium citrate antigen	Servicebio	G1202
retrieval solution (pH 6.0)		
PBS solution	Servicebio	G0002
3%H ₂ O ₂	SCRC	
FITC-TSA	Servicebio	G1222
CY3-TSA	Servicebio	G1223
647-TSA	Servicebio	G1224
Spontaneous fluorescence	Servicebio	G1221
quenching reagent		
BSA	Servicebio	G5001
Primary antibody		
Secondary antibody		
DAPI	Servicebio	G1012

2. Procedure

2.1 Deparaffinize and rehydrate: Incubate sections in 3 changes of Biodewax and Clear Solution, 10 min each. Dehydrate in 3 changes of pure ethanol for 5 min. Wash in distilled water.

2.2 Antigen retrieval: Immerse the slides in EDTA antigen retrieval buffer (pH 8.0) and maintain at a sub-boiling temperature for 8 min, standing for 8 min and then followed by another sub-boiling temperature for 7 min. Be sure to prevent buffer solution evaporate. Let air cooling. Wash three times with PBS (pH 7.4) in a Rocker device, 5 min each. Use the right antigen retrieval buffer and heat extent according to tissue characteristics.

2.3 Circle, Block endogenous peroxidase: Wash three times with PBS (pH 7.4) in a Rocker device, 5 min each eliminate obvious liquid, mark the objective tissue with liquid blocker pen. Immerse in 3% H₂O₂ and incubate at room temperature for 15 min, kept in dark place. Then wash again three times with PBS (pH 7.4) in a Rocker device, 5 min each.

2.4 Block with serum: Eliminate obvious liquid, mark the objective tissue with liquid blocker pen. Cover objective tissues with 10% donkey serum (for the case of primary antibody originated from goat) or 3% BSA (for the case of primary antibody originated from others) at room temperature for 30 min.

2.5 First primary antibody: Throw away the blocking solution slightly. Incubate slides with the first primary antibody (diluted with PBS appropriately) overnight at 4 °C, placed in a wet box containing a little water.

2.6 Corresponding secondary antibody marked with HRP: Wash slides three times with PBS (pH 7.4) in a Rocker device, 5 min each. Then throw away liquid slightly. Cover objective tissue with secondary antibody (appropriately respond to first primary antibody in species), incubate at room temperature for 50 min in dark condition.

2.7 CY3-TSA solution: Wash slides three times with PBS (pH 7.4) in a Rocker device, 5 min each. Incubate slides with CY3-TSA solution (diluted with TBST appropriately) for 10 min in dark condition. After that, wash slides three times with TBST in a Rocker device, 5 min each.

2.8 Microwave treatment: Immerse the slides in EDTA antigen retrieval buffer (pH 8.0) and maintain at a sub-boiling temperature for 8 min, standing for 8 min and then followed by another sub-boiling temperature for 7 min, to remove the primary antibodies and Secondary antibodies combined with tissue. Be sure to prevent buffer solution evaporate.

2.9 Second primary antibody: Incubate slides with primary antibody (diluted with PBS appropriately) overnight at 4 °C, placed in a wet box containing a little water.

2.10 Second corresponding secondary antibody marked with HRP: Wash slides three times with PBS (pH 7.4) in a Rocker device, 5 min each. Then throw away liquid slightly. Cover objective tissue with secondary antibody (appropriately respond to second primary antibody in species), incubate at room temperature for 50 min in dark condition.

2.11 FITC-TSA solution: Wash slides three times with PBS (pH 7.4) in a Rocker device, 5 min each. Incubate slides with FITC-TSA solution (diluted with TBST appropriately) for 10 min in dark condition. After that, wash slides three times with TBST in a Rocker device, 5 min each.

2.12 Microwave treatment: Immerse the slides in EDTA antigen retrieval buffer (pH 8.0) and maintain at a sub-boiling temperature for 8 min, standing for 8 min and then followed by another sub-boiling temperature for 7 min, to remove the primary antibodies and Secondary antibodies combined with tissue. Be sure to prevent buffer solution evaporate.

2.13 Incubate the third primary antibody: Incubate slides with primary antibody (diluted with PBS appropriately) overnight at 4 °C, placed in a wet box containing a little water.

2.14 The third corresponding secondary antibody marked with HRP: Wash slides three times with PBS (pH 7.4) in a Rocker device, 5 min each. Then throw away liquid slightly. Cover objective tissue with secondary antibody (appropriately respond to second primary antibody in species), incubate at room temperature for 50 min in dark condition.

2.15 647-TSA solution: Wash slides three times with PBS (pH 7.4) in a Rocker device, 5 min each. Incubate slides with 647-TSA solution (diluted with TBST appropriately) for 10 min in dark condition. After that, wash slides three times with TBST in a Rocker device, 5 min each.

2.16 Microwave treatment: Immerse the slides in EDTA antigen retrieval buffer (pH 8.0) and maintain at a sub-boiling temperature for 8 min, standing for 8 min and then followed by another sub-boiling temperature for 7 min, to remove the primary antibodies and Secondary antibodies

combined with tissue. Be sure to prevent buffer solution evaporate.

2.17 Incubate the fourth primary antibody: Incubate slides with primary antibody (diluted with PBS appropriately) overnight at 4 °C, placed in a wet box containing a little water.

2.18 The fourth corresponding secondary antibody marked with 594 labeled Fluorescence secondary antibody: Wash slides three times with PBS (pH 7.4) in a Rocker device, 5 min each. Then throw away liquid slightly. Cover objective tissue with secondary antibody (appropriately respond to second primary antibody in species), incubate at room temperature for 50 min in dark condition.

2.19 DAPI counterstain in nucleus: Incubate with DAPI solution at room temperature for 10 min, kept in dark place.

2.20 Spontaneous fluorescence quenching: Eliminate obvious liquid, incubate slides with spontaneous fluorescence quenching reagent for 5 min, then wash slides under flowing water for 10 min.

2.21 Mount: Wash three times with PBS (pH 7.4) in a Rocker device, 5 min each. Throw away liquid slightly, then coverslip with anti-fade mounting medium.

2.22 Detection and collect images by slice scanner: DAPI glows blue by UV excitation wavelength 330-380 nm and emission wavelength 420 nm; FITC glows green by excitation wavelength 465-495 nm and emission wavelength 515-555 nm; CY3 glows red by excitation wavelength 510-560 nm and emission wavelength 590 nm. 647 glows pink by excitation wavelength 608-648 nm and emission wavelength 672-712 nm. (647 was originally red, in order to distinguish it from CY3, we set it to pink light.). 594 glows purplish red by excitation wavelength 594 nm and emission wavelength 615 nm. (594 was originally red, in order to distinguish it from CY3, we set it to purplish red light.)

3. Introduction to Technical Principles

The main principle is based on Tyramide Signal Amplification (TSA), hereinafter referred to as TSA technology. TSA is an enzyme detection method based on the catalytic activity of horseradish peroxidase HRP for high-density in situ labeling of target proteins or nucleic acids. The main principle of the technology is that fluorescent-labeled tyramine becomes activated tyramine under the action of HRP and H₂O₂, which attaches to the protein tyrosine residues around the target. This binding is covalent binding, while the primary antibody is non-covalent binding with the target, and the secondary antibody is non-covalent binding with the primary antibody. Through microwave repair treatment, the first round of primary antibody and secondary antibody are eluted, and fluorescently labeled tyramine is still attached around the target. When detecting the second target, it is equivalent to a completely new round of markers, without considering whether the antibody of the second round cross-reacts with the antibody of the first round. Multiple targets can

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Wuhan Servicebio Technology Co., Ltd.

be labeled only by changing different fluorescence markers. Double or multiple fluorescent staining was performed using different fluorescent tyramines by repeated immunolabeling.

4. Results

Nucleus is blue by labeling with DAPI. Positive cells are green, red, pink or purplish red according to the fluorescent labels used.

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