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Immunofluorescence Protocol (Frozen-slides)

1. Apparatus and Reagents

1.1 Apparatus

Name	Producer	Model
Freezing microtome	Thermo	Cryotome E
Glass microscope slides	Servicebio	
Coverslips	CITOTEST	10212432C
Microwave	Glanze	P70D20TL-P4
Rocker	Servicebio	TSY-B
Vortex	Servicebio	MX-F
Micro-centrifuge	Servicebio	D1008E
Pipettor	Dragon	KE0003087/KA0056573
Liquid blocker pen	Servicebio	G6100
Refrigerator	Haier	BCD-192TGN
Ortho-Fluorescent Microscopy	Nikon	NIKON ECLIPSE C1
Imaging system	Nikon	NIKON DS-U3

1.2 Major reagents

Name	Producer	Code	Dilution
OCT embedding medium	Servicebio	G6059-110ML	
EDTA antigen retrieval solution (pH 8.0)	Servicebio	G1206	
EDTA antigen retrieval solution (pH 9.0)	Servicebio	G1203	
Sodium citrate antigen retrieval solution (pH 6.0)	Servicebio	G1202	
BSA	Servicebio	G5001	
4% paraformaldehyde	Servicebio	G1101	
PBS solution	Servicebio	G0002	
Spontaneous fluorescence quenching reagent	Servicebio	G1221	

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Add: 5th Floor, C22 Building, Biopark, No. 388 Gaoxin 2nd Road, East Lake High-tech Developing Zone,
Wuhan, Hubei, China, 430079 Website: <http://www.servicebio.com> Tel: 027-5111-3188

Primary antibody

Fluorescent-labelled

secondary antibody

DAPI	Servicebio	G1012
anti-fade mounting medium	Servicebio	G1401

2. Procedure

2.1 Fix frozen-slides: Frozen slides are baked in a 37°C oven for 10-20 minutes to control the moisture. Fix in paraformaldehyde 30 min and then dry in air. Wash three times with PBS (pH 7.4) in a Rocker device, 5 min each.

2.2 Antigen retrieval: Immerse the slides in EDTA antigen retrieval buffer (pH 8.0) and maintain at a sub-boiling temperature for 8 min, standing for 8 min and then followed by another sub-boiling temperature for 7 min. Be sure to prevent buffer solution evaporate. Let air cooling. Wash three times with PBS (pH 7.4) in a Rocker device, 5 min each. Use the right antigen retrieval buffer and heat extent according to tissue characteristics.

2.3 Circle and Serum blocking: Eliminate obvious liquid, mark the objective tissue with liquid blocker pen. Add 3% BSA to cover the marked tissue to block non-specific binding for 30 min. Cover objective area with 10% donkey serum (for the case of primary antibody originated from goat) or 3% BSA (for the case of primary antibody originated from others)

2.4 Primary antibody: Throw away the blocking solution slightly. Incubate slides with primary antibody (diluted with PBS appropriately) overnight at 4 °C, placed in a wet box containing a little water.

2.5 Secondary antibody: Wash slides three times with PBS (pH 7.4) in a Rocker device, 5 min each. Then throw away liquid slightly. Cover objective tissue with secondary antibody (appropriately respond to primary antibody in species), incubate at room temperature for 50 min in dark condition.

2.6 DAPI counterstain in nucleus: Wash three times with PBS (pH 7.4) in a Rocker device, 5 min each. Then incubate with DAPI solution at room temperature for 10 min, kept in dark place.

2.7 Spontaneous fluorescence quenching: Wash three times with PBS (pH 7.4) in a Rocker device, 5 min each. Add spontaneous fluorescence quenching reagent to incubate for 5 min. Wash in running tap water for 10 min.

2.8 Mount: Throw away liquid slightly, then coverslip with anti-fade mounting medium.

2.9 Microscopy detection and collect images by Fluorescent Microscopy. DAPI glows blue by UV excitation wavelength 330-380 nm and emission wavelength 420 nm; FITC glows green by excitation wavelength 465-495 nm and emission wavelength 515-555 nm; CY3 glows red by

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excitation wavelength 510-560 nm and emission wavelength 590 nm.

3. Results

Nucleus is blue by labeling with DAPI. Positive cells are green or red according to the fluorescent labels used.

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Wuhan Servicebio Technology Co., Ltd.