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Add: 5th Floor, C22 Building, Biopark, No. 388 Gaoxin 2nd Road, East Lake High-tech Developing Zone, Wuhan, Hubei, China, 430079 Website: <http://www.servicebio.com> Tel: 027-5111-3188

蛋白酶 K	Servicebio	G1234
DAPI	Servicebio	G1012
抗荧光淬灭封片剂	Servicebio	G1401
杂交缓冲液	Servicebio	G3016-3
anti-DIG-HRP	Jackson	200-032-156
FITC-Tyramide	Servicebio	G1222

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1. Apparatus and reagents

1.1 Apparatus

Name	Producer	Model
Dehydrator	DIAPATH	Donatello
Glass slide	Servicebio	
Coverslip	Citotest	10212432C
Pathologic microtome	Leica	RM2016
Enzyme-free centrifuge tubes	Servicebio	EP-150-M
Shaker	Servicebio	TSY-B
Vortex	Servicebio	MV-100
Pipettor	Dragon	KE0003087/KA0056573
Liquid Blocker PAP Pen	Servicebio	G6100
Microscopy	Nikon	Nikon Eclipse ci
Imaging system	Nikon	Nikon DS-U3
Incubator	LABOTERY	GSP-70
autoclave	Panasonic	MLS-3751L-PC
Paraffin embedding machine	WHJJ	JB-P5
Frozen flat	WHJJ	JB-L5
Tissue spreader	Kedee	KD-P

1.2 Major reagents

reagent	manufacturer	article number
DEPC	Amresco	E174
4% of paraformaldehyde (DEPC water)	Servicebio	G1113
Paraffin wax	Sakura	G1113
Ethanol	SCRC	100092683
BioDewax and Clear solution	Servicebio	G1128
PBS	Servicebio	G0002
20×SSC solution	Servicebio	G3016-4
BSA	Servicebio	G5001

Proteinase K	Servicebio	G1234
DAPI	Servicebio	G1012
Anti-fluorescence quenching sealing tablets	Servicebio	G1401
hybridization buffer	Servicebio	G3016-3
anti-DIG-HRP	Jackson	G3016-3
FITC- Tyramide	Servicebio	G1222

2. The steps of the experiment

2.1 Organization fixed: take out the organization, wash clean, then Immediately put in the fixed fluid (DEPC) above 12h.

2.2 Dehydration: The tissue is dehydrated by gradient alcohol, paraffin, embedding.

2.3 Section: the paraffin is sliced through the slicer, the piece of the slice machine and the 62 - degree oven roast for 2 hours.

2.4 Dewaxing and dehydration: soak sections in 2 changes of Dewaxing Transparent Liquid, 15 minutes each. Dehydrate in 2 changes of pure ethanol for 5 minutes each. Then, followed respectively by dehydrating in gradient ethanol of 85%and 75% ethanol 5 minutes for each. Wash in DEPC dilution.

2.5 Digestion: according to the tissue fixation time, the slices are boiled in the retrieval solution for 10-15 minutes and naturally cooled. Mark the objective tissue with liquid blocker pen, according to the characteristics of tissues, Add proteinase K(20 ug/ml) working solution to cover objectives and incubate at 37°C for___min. Washing with pure water, then wash three times with PBS (pH 7.4) in a Rocker device, 5 min each.

2.6 Block endogenous peroxidase: add 3% methanol-H₂O₂, incubate in dark at room temperature for 15min. Wash slides in PBS(PH7.4) three times for 5 min each with gentle agitation.

2.7 Pre-hybridization: add Pre-hybridization solution to each section and incubate for 1 h at 37°C.

2.8 Hybridization: remove the pre-hybridization solution, add the probe hybridization solution with concentration of ___, and incubate the section in a humidity chamber and hybridize overnight at___°C.

2.9 Washing: remove the hybridization solution. Wash sections with 2×SSC for 10 min at 37°C , Wash sections in 1×SSC two times for 5 min each at 37°C, and wash in 0.5×SSC for 10 min at room temperature. Formamide washing can be added if there are more non-specific hybrids.

2.10.Imaging oligo (DIG): Add the hybridization solution containing Imaging Oligo (DIG), Dilution ratio 1:400. 42 ° FOR 3 hours. Wash sections with 2×SSC for 5 min at 37°C , Wash sections in 1×SSC two times for 5 min each at 37°C, and wash in 0.5×SSC for 10 min at room

temperature.

2.11 **Blocking:** add blocking solution (rabbit serum) to the section and incubate at room temperature for 30 min.

2.12 Add the mouse anti-digoxigenin-labeled peroxidase (anti-DIG-HRP): remove the blocking solution and add anti-DIG-HRP. Incubate at 37 °C for 40 min, then Wash sections in PBS four times for 5 min each.

2.13 Tyramide **developing:** dry sections slightly, add fresh prepared Tyramide chromogenic reagent to marked tissue. Reaction in dark for 5 min at room temperature. Then wash sections in PBS three times for 10 min each.

2.14 **Stain cell nuclei (counter stain):** incubate with DAPI for 8min in the dark, and then mounting.

2.15 **Microscopic examination and photography:** to take photos with positive fluorescence microscope. DAPI glows blue by UV excitation wavelength 330-380 nm and emission wavelength 420 nm; FAM glows green by excitation wavelength 465-495 nm and emission wavelength 515-555 nm; CY3 glows red by excitation wavelength 510-560 nm and emission wavelength 590 nm.

3. Interpretation of the results

The nuclear stained by DAPI were blue under ultraviolet excitation, and the positive expression was a kind of fluorescence labeled by corresponding luciferin. FAM (488) appears green, cy3 appears red. The results of mRNA in situ hybridization were cytoplasmic positive and a few nuclear positive were normal. MicRNA and lncRNA were expressed differently. According to the expression, Different fluorescence brightness is strong or weak.

Note: All reagents, instrument need RNase free.

Attached table 1 probe information.

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Wuhan Servicebio Technology Co., Ltd.